WELCOME TO ALPS® SOUTH

ALPS® South is an American company which has been specializing in the production of prosthetic liners since 1993. Since then we have expanded globally, with offices in 5 countries, a consolidated network of agents and over 250 trusted distributors spread out all over the world.

Thanks to our dedication to innovation, coupled with our eagerness to interpret and even anticipate the needs of prosthetic technicians and wearers, ALPS® has gone far beyond our initial production only of prosthetic liners, and developed a whole range of high-quality prosthetic solutions, including more than 75 registered patents spanning an array of product categories. This catalogue makes it easy to access information about the entire range of ALPS® products to identify the best prosthetic solutions based on both the lifestyle and clinical profile of each individual patient.

PREMIUM QUALITY | 100% MADE IN USA

Investment in innovation, research and quality, both in terms of products and working methods, continues to be a priority here at ALPS®, allowing us to offer premium solutions. The entire production process, from the knitting of fabrics to the manufacturing of the liners, takes place at our St. Petersburg, Florida plant. This means we can assure maximum quality control at every step of production. In addition ALPS® guarantees the quality and safety of its products with certification from some of the most important internationally recognized authorizing institutions. ALPS® has been manufacturing FDA-approved medical devices since 1989. ALPS® products are designed and produced according to ISO13485 and GMP (Good manufacturing practice) standards, meeting the requirements of the EU regulations CEE 93/42 UE and conforming to MDR 2017/745.





ALPS® SOUTH MAKING LIVES BETTER™

"Making Lives Better™" is the mission which has been driving the company since 1988, inspiring our job every day and pushing us to provide superior products for you and for every prosthetics user.

We pride ourselves on our fair business practices with customers, vendors, and employees alike, while maintaining a brand image that makes our staff, collaborators and end-users of our products all proud to be associated with ALPS®.

LINER SPECIALISTS SINCE 1993

When it comes to optimizing prosthetic fit and patient outcomes, liner selection plays a critical role. With a combination of five elastomers, five fabrics, and five suspension options, ALPS offers one of the most versatile and high-performing liner portfolios on the market, designed to meet the diverse needs of patients and technicians.

Throughout the catalog, you'll find a QR code that links to a simple, easy-to-understand overview of the prosthetic liner technologies we offer. This webpage will offer helpful explanations of each elastomer, fabric, and suspension type, empowering you to make informed decisions and select the most appropriate liner solution for every patient.

Liner Technology

5 Elastomers. 5 Fabrics. 5 Suspensions.



INDEX CATALOG

LINERS OVERVIEW	8
SLEEVES OVERVIEW	14
HOW TO SIZE LINERS AND SLEEVES	
ALPS TECHNOLOGIES	18
GEL TECHNOLOGY	20
KNITTING TECHNOLOGY	22
SILICONE TECHNOLOGY	24
NRG TECHNOLOGY®	26
DISTAL END CONSTRUCTION	28
HOW ANTIOXIDANTS WORK	30
CASE STUDIES	
PROSTHETIC LINERS	34
ProFit Liner	36
SoftSil® Liner	38
UltraSeal Silicone Liner	40
UltraSeal Gel Liner	42
EasyLiner®	44
Silicone PRO Liner	46
Superior Performance Liner	48
Superior Performance HD Liner	50
General Purpose Liner	52
Eco™ Liner	54
Smart Seal™ Liner	56
Extreme™ Liner	58
Extreme TM HD Liner	14



EasyFlex™ Liner	62
EasyLiner® SuperStretch	64
Anterior Posterior Tapered Liner	66
PSI Liner	68
Liberty™ Liner	70
PROSTHETIC SUSPENSION SYSTEMS	72
SLEEVES	
Reinforced Flex Sleeve	74
Reinforced Sleeve	
Flex Sleeve	78
EasySleeve	80
EasySleeve SuperStretch	82
Superior Suspension Sleeve	84
Superior Performance Sleeve	86
Extreme Sleeve	88
V.I.V.A. Sleeve	90
Anterior Posterior Tapered Sleeve	
EasySleeve Fabric Reinforced	94
Vacuum Integrated Pump	96
LOCKS	
VIP® Lock	98
Lock with Pyramid	100
Lock with Flat Bottom	102
VOLUME CHANGE MANAGEMENT	
Prosthetic Socks	106
Skin Reliever	108



INDEX CATALOG

LINER AND RESIDUAL LIMB CARE

Antioxidant Cream	110
Silicone Skin Lotion	11
Prosthetic Lubricant	
Prosthetic Cleanser	
Prosthetic Anti-Perspirant	
Prosthetic Ointment	
EXTRA	
PadSil	
PROSTHETIC FEET	118
BioStep® EVO	120
BioStep® PRO	122
BioStep® CL	124
BioStep® XT	126
BioStep® LP	128
PROSTHETIC KNEE	
ProStride™	130
ALPS® WARRANTY POLICY	132
CONTACTS	13/





Locking Liner



Cushion Liner





EasyGel: contains antioxidants that are an effective scavenger of free radicals and hydroxyl groups. Best for use on amputees with adherent scar tissue, skin irritation, sensitive skin and diabetics.



Grip Gel: tackiness helps in reducing bunching in the popliteal area. Best for use on amputees with healthy skin tissue.



High Density Gel: the relatively firm nature of the HD Gel combines comfort with control and stability. Best for use on patients who are transitioning from a silicone liner to a gel liner.



OptiGel: ALPS® exclusive, soft, white gel that ensures excellent adaptation to the residual limb shape, making it the optimal choice for users who demand ultimate comfort and control.



Silicone



OptiSil (White Silicone)



In-house knitted fabric with high quality yarns



Black Fabric with highly reduced vertical stretch



Brown Fabric with reduced vertical stretch



Beige Fabric with both horizontal and vertical stretch



Gray Fabric with increased reduction of vertical elongation

CTIVITY LEVE



K1: low activity



K3: middle/high activity



K2: middle/low activity



K4: high activity



LINERS OVERVIEW

PICTURE	NAME	PART #	SUSPENSION	MATERIAL
a profit ◆	ProFitLiner	PFGDT PFGFR	UŲ	OptiGel
1 Xerra	SoftSil® Liner	SSGSL SSGSC	UU	OptiSil
1 American	UltraSeal Silicone Liner	USHSL USHSC USLSL USLSC USUSL USUSC	UU	OptiSil
1 Accessor	UltraSeal Gel Liner	USHGL USHGC USLGL USLGC USUGL USUGC	UŲ	OptiGel
	EasyLiner®	ELDT ELFR	UU	EasyGel
	Silicone Pro Liner	LSC LSL	UŲ	Silicone
0.7 	Superior Performance Liner	SPDT SPFR	U	GripGel



THICKNESS	ACTIVITY LEVEL	RECOMMENDED FOR	PAGE
3 mm 6 mm		Patients who demand comfort, control, and stability needed for everyday use or high performance activities.	36
2 mm	\$\frac{\cdots}{\cdots} \frac{\cdots}{\cdots} \frac{\cdots}{\cdots	Patients with medium to high activity levels as the fabric offers a high degree of control, while the soft silicone provides comfort. Further recommended for patients living in hot areas.	38
2 mm	\(\hat{\chi}\) \(\hat{\chi}\) \(\hat{\chi}\) \(\hat{\chi}\) \(\hat{\chi}\) \(\hat{\chi}\) \(\kappa\) \(\kappa	Transtibial and transfemoral users with medium to long residual limbs and for knee disarticulation.	40
3 mm	パ	Transtibial and transfemoral users with medium to long residual limbs and for knee disarticulation.	42
3 mm 6 mm 6/3 mm distal proximal taper	\(\hat{\chi}\) \(\hat{\chi}\) \(\hat{\chi}\) \(\hat{\chi}\) \(\hat{\chi}\) \(\hat{\chi}\) \(\kappa\) \(\kappa	Diabetic patients and those with sensitive skin and whenever the comfort and the protection of the residual limb is a priority.	44
2 mm	当 学 学 K1 × K2 × K3 × K4 ×	Residual limbs with redundant soft tissue and for transfemoral amputees. Further recommended to counter perspiration issues and for patients living in hot areas.	46
3 mm 6 mm	Å Å Å Å Å Å Å Å Å Å	Transfemoral and transtibial patients with moderate activity levels, and good skin condition who are looking for an improved control of the prosthetic leg.	48



LINERS OVERVIEW

PICTURE	NAME	PART #	SUSPENSION	MATERIAL
	Superior Peformance HD Liner	SPDTHD SPFRHD	UU	High Density Gel
	General Purpose Liner	GPDT GPFR	UU	GripGel
	Eco™ Liner	ECDTHD ECFRHD	UŲ	High Density Gel
	Smart Seal™ Liner	SADT SAFR SBFR	UŲ	High Density Gel
	Extreme™ Liner	AKDT AKFR	UU	GripGel
	Extreme™ HD Liner	AKDTHD AKFRHD	UU	High Density Gel
	EasyFlex™ Liner	EZF	Ų	EasyGel



THICKNESS	ACTIVITY LEVEL	RECOMMENDED FOR	PAGE
3 mm 6 mm	\$\frac{\cdots}{\cdots} \frac{\cdots}{\cdots} \frac{\cdots}{\cdots	Patients with medium and high activity levels who require both comfort and elevated levels of control of the prosthetic leg. It is also recommended for patients whose residual limbs have excessive soft tissue and for transfemoral amputees.	50
3 mm 6 mm	パ パ ペ ペ K1 × K2 × K3 × K4 ×	Patients with low activity levels and good skin conditions; it is also indicated for irregularly shaped residual limbs and those with bony protuberances.	52
3 mm 6 mm	労 	Patients with medium to high activity levels, whether with transtibial or transfemoral prostheses, and for those whose residual limbs present redundant soft tissues.	54
3 mm 6 mm	外	Transtibial and transfemoral users with medium to long residual limbs and for knee disarticulation.	56
3 mm 6 mm	Å Å Å Å Å K 1 ∨ K 2 ∨ K 3 ∨ K 4	Patients, both transfemoral and transtibial, with medium activity levels whose skin is in good condition and who are aiming to maintain a good level of control of the prosthetic leg.	58
3 mm 6 mm	Å Å Å Å Å K 4 ✓ K 3 ✓ K 4 ✓	Patients with medium and high activity levels as it offers a high level of control while preserving comfort. It is also recommended for patients whose residual limbs have excessive soft tissue and for transfemoral amputees.	60
3 mm 6 mm	パ パ ペ ペ K1 × K2 × K3 × K4 ×	It is recommended for diabetic and sensitive skin transtibial patients who have medium to high activity levels and using a pin and lock suspension.	62



LINERS OVERVIEW

PICTURE	NAME	PART#	SUSPENSION	MATERIAL
	EasyLiner® SuperStretch	ELPX	U	EasyGel
	Anterior Posterior Tapered Liner	APDT APFR	U	High Density Gel
÷	PSI Liner	PSI	U	EasyGel
4	Liberty™ Liner	VSDT	Ų	High Density Gel



THICKNESS	ACTIVITY LEVEL	RECOMMENDED FOR	PAGE
6 mm	№ № № № № № №	It is recommended for diabetic patients and those with sensitive skin or irregularly shaped residual limbs.	64
6/3 mm tapered	労	It is indicated for transtibial users with medium and high activity levels who desire greater control and freedom of movement without sacrificing comfort.	66
6 mm 6/3 mm distal proximal taper	№ № № № № № №	It is recommended for transtibial patients with low and medium activity levels and sensitive skin.	68
3 mm 6 mm	\$\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\	It is recommended to counter the problem of excessive perspiration.	70



SLEEVES OVERVIEW

PICTURE	NAME	PART #	MATERIAL	THICKNESS	KEY FEATURES
AX:	Reinforced Flex Sleeve	RFX®	OptiGel	3 mm	30° pre-flexion + an added layer of internal fabric for increased durability
₹FS man	Reinforced Sleeve	RFS	OptiGel	3 mm	Constructed with an extra layer of internal fabric to deliver superior durability
	Flex Sleeve	SFX	HD Gel	3 mm	30° pre-flexed
	Easy Sleeve	SFR	EasyGel	3 mm 6 mm	For sensitive skins
	EasySleeve SuperStretch	SLE	EasyGel	6 mm	Without fabric for freedom of movement
	Superior Suspension Sleeve	SFB SFBHD	HD Gel GripGel	2 mm	2 mm thinner profile for less bulk
	Superior Performance Sleeve	SFP SFPHD	HD Gel GripGel	3 mm 6 mm	The innovative black fabric provides comfort and stability



ACTIVITY LEVEL

















HOW TO SIZE LINERS



Measure the circumference of the residual limb at 6 cm above the distal end. Use the chart to select the size of the liner.

Example: for a measured circumference of 30 cm, select a 28 size liner.

HOW TO SIZE SLEEVES



Measure the circumference of the leg at 10 cm above the center of the knee. Use the chart to select the size of the sleeve.

Example: for a measured circumference of 32 cm, select a 24 size sleeve.



SLEEVES OVERVIEW

PICTURE	NAME	PART #	MATERIAL	THICKNESS	KEY FEATURES
	Extreme Sleeve	SFK	GripGel	3 mm 6 mm	Limited vertical stretch, ideal for active patients
	V.I.V.A	VIVA	GripGel	3 mm 6 mm	The built-in valve evacuates the air in the socket
	Anterior Posterior Tapered Sleeve	APS	HD Gel	6/3 mm tapered	Gradually tapered from the anterior to the posterior to provide superior comfort
	EasySleeve Fabric Reinforced	SFS	GripGel	3 mm 6 mm	Achieves secure suspension by sealing with the skin without restricting circulation



ACTIVITY LEVEL






K1 ✓ K2 ✓ K3 K4






HOW TO SIZE LINERS



Measure the circumference of the residual limb at 6 cm above the distal end. Use the chart to select the size of the liner.

Example: for a measured circumference of 30 cm, select a 28 size liner.


HOW TO SIZE SLEEVES




Measure the circumference of the leg at 10 cm above the center of the knee. Use the chart to select the size of the sleeve.


Example: for a measured circumference of 32 cm, select a 24 size sleeve.








GEL TECHNOLOGY



ALPS® expertise in the field of chemistry and our position as a leading manufacturer of advanced gel-based medical devices have led to the development of a wide range of prosthetic liners and sleeves over the years, characterized by an innovative and unique product: our patented gel.

ALPS® Gel offers an optimum level of protection against shear forces and other stresses, guaranteeing the perfect combination of control and uniform distribution of pressure within the socket thanks to its characteristic low modulus of elasticity.

ALPS® Gel is far more flexible and resistant with respect to other thermoplastic gels being used in the lower limb prosthetic liner industry. ALPS® Gel offers up to six times the elasticity of other thermoplastic (TPE) gels on the market. It is this important mechanical property that allows ALPS® Gel to guarantee prosthetic leg wearers the comfort, durability, control and safety they need to grant them the mobility and freedom of movement.







Four types of ALPS® Gels to serve the mobility needs of prosthetic leg wearers.

ALPS® makes four types of gels to meet the needs of every individual prosthetic leg user whose requirements vary based on their clinical condition and lifestyle: EasyGel, HD Gel, GripGel, and OptiGel. The full range of ALPS® Gels means that each patient will be able to find his or her own, unique, perfect comfort.



EASYGEL

EasyGel was the first kind of gel developed by ALPS®. It is characterized by its soft, powdery texture. It is rich in antioxidants that are an effective scavenger of free radicals and hydroxyl groups. Best for use on amputees with adherent scar tissue, skin irritations, sensitive skin and diabetics.



HIGH DENSITY GEL (Silicone equivalent)

Dense, firm texture. Relatively firm nature of the HD Gel combines comfort with control and stability. Best for use on highly active patients who are transitioning from a silicone liner to a gel liner.



GRIPGEL

Smooth, tacky texture. Tackiness helps in reducing bunching in the popliteal area. Best for use on amputees with healthy skin tissue.



OPTIGEL

This exlusive, white gel ensures excellent adhesion to the residual limb, making it the optimal choice for users who demand ultimate comfort, control, and stability.



KNITTING TECHNOLOGY

The choice of fabric used for producing prosthetic sleeves and liners is fundamental for guaranteeing maximum comfort and safety.

This is why ALPS® began manufacturing its own fabrics years ago: to be sure to offer a product that best meets the functional needs of prosthetic wearers. The fabric needs to do more than just protect the materials used in making the liners and extend their durability; it must also be able to generate a greater, lesser or differentiated level of tension depending on the needs of the individual user. ALPS® manufactures and uses five types of fabric: beige, brown, black, gray and knitted.

These five types of fabric are not differentiated merely by color and manufacturing process, they also each have distinctive characteristics and functions in order to provide each prosthetic leg wearer with the specific benefits he or she needs based on individual lifestyles and clinical conditions.







BEIGE FABRIC

Beige fabric stretches both vertically and horizontally

The fabric's stretch allows the liner to adapt to the shape of the residual limb as it moves to guarantee maximum comfort.



BROWN FABRIC

Brown fabric has reduced vertical elasticity

Vertical stretch is reduced by 80% compared to the beige fabric, allowing a greater level of control and stability of the prosthesis. Ideal for active patients.



BLACK FABRIC

Black fabric with greatly reduced vertical stretch

Vertical stretch is 90% lower than the beige fabric, allowing a greater level of control and stability of the prosthesis. Ideal for active patients.



GRAY FABRIC

Gray fabric with increased reduction of vertical stretch

Minimal vertical stretch allows for a high level of control and freedom of movement. Ideal for active patients.



KNITTED FABRIC

ALPS® years of experience in manufacturing fabrics has played a fundamental role in the development of our latest innovation: products created using ALPS® knitted fabric.

Our latest fabric manufacturing process has allowed us to:

- Create products totally devoid of any stitching
- Maintain control of the tension in each section of the liner to maximize and reduce the effects of pistoning
- Integrate a matrix in the distal construction of our locking liners
- Integrate a distal matrix to limit elongation / U matrix or Regular Matrix



SILICONE TECHNOLOGY

ALPS® has introduced innovative and skin-friendly silicone prosthetic liners for lower limb amputation.


These liners provide a solution to some of the clinical problems and disadvantages of silicone in skin-contact applications due to the inherently high coefficient of static friction typical of silicone polymers.

ALPS® expertise in processing silicone allowed our engineers to come up with a new manufacturing method which improved the friction characteristics of the silicone prosthetic liner.

The surface of the prosthetic liner has been modified, to create micro craters that reduce the static friction without any additives, resulting in a completely skin-friendly silicone.

Silicone with the modified surface layer shows a reduction of 80% of the coefficient of static friction, as shown in the graph below.

Friction Testing







ALPS® Silicone liners are manufactured using pure silicone. This makes our silicone liners suitable for a wide range of amputees, including elderly and vascular patients with a lower limb amputation.

The silicone products stand alongside our other liners manufactured in EasyGel, Grip Gel, High Density Gel (Silicone Equivalent), and the newest OptiGel.

ALPS®' silicone is lightweight, hypoallergenic, cool and additive-free and creates only one-fifth the friction of normal silicone. The technology used to create this particular silicone is of fundamental importance, allowing ALPS® to design and produce silicone liners for prosthetic legs which can guarantee:

- Freedom of movement
- Maximum gentleness on the skin
- Cooling effect
- · Ease of cleaning





NRG TECHNOLOGY®



Inspired by Nature.

Empowered by Technology.


NRG - Natural Reactive Gait - is a patented innovative dynamic response foot technology designed to replicate the natural movement of the calf and shin muscles.





ONE FOOT, MULTIPLE POSSIBILITIES.





KEY FEATURES



Unique multilayer construction



Exclusive and patented reverse shank design



Advanced flexible ergonomic design

KEY BENEFITS

- Physiological and energy efficient gait
- Smooth rollover during the whole gait cycle
- High dynamic response to support an active lifestyle
- Possibility to do a wide range of activities with just one single foot
- Ground compliance and stability on all terrains
- Reduction of secondary issues thanks to a more balanced posture and symmetrical gait
- Low maintenance
- Lightweight
- Waterproof






DISTAL END TECHNOLOGY

Our distal end construction delivers a superior, high-quality product that provides the comfort, protection, and durability essential for a successful prosthetic liner.

How we do it.

Our team of engineers works meticulously to ensure that every distal end is constructed to the highest standards, using only the best materials and components.

That's why our prosthetic liners are designed using an **internal distal end construction method** which fuses the three primary liner materials- gel, fabric, and umbrella- on the inside of the liner, instead of being glued on from the outside like many liners available today.



Integrating these materials into a cohesive unit creates a secure anchoring point for the umbrella, ensuring that the liner can easily connect to a prosthetic device without fear of damage or detachment. Not only does this guarantee the reliability of the entire system, but it also provides unparalleled comfort and performance.

4 distal end styles and 6 different construction processes.



Each prosthetic liner falls under 1 of 4 distal end styles, which is determined by umbrella type. The liners within each style are then further grouped based on its internal construction process (i.e., the components used during assembly).

Each of the six construction processes is unique and complex, involving several elements that merge together to create a highly durable structure that caters to the liner's purpose and recommended use.


Refer to the Table on page 29 to view the specs of each distal end style and the components used in their construction.



Distal E	nd Styles	Distal End	Construction		Distal End Co	mponents	
Style Name	Umbrella Type	Internal Construction	Liner(s)/Construction Process	Grommet	Internal Reinforcement Fabric	U Band	Encapsulation
Easy	Wide, Flat	/	EasyLiner	✓	✓		Plastic
Retro	Standard	,	GP, AK, EZF, VSDT, SP	/	✓		Silicone
Ketro	Standard	✓	Eco	/		/	Silicone
01:	Low-Profile	,	ProFit	/		/	Aluminum*
Choice	Low-Profile	~	SoftSil			/	Aluminum*
			UltraSeal			/	Aluminum*
Classic	Low-Profile	/	Silicone Pro			/	

Gel and fabric type not included as it may vary.

*U.S. Patent No. 1,076,089



Deep Dive

Style: Choice Construction Process: Pro-Fit

A. Gel

The thickness layer of the OptiGel increases at the distal end to provide cushion. It also embeds the base of the umbrella (B) and helps secure the fabric layers (D, E) in place.

B. Umbrella

An internal metal umbrella makes the liner more secure and durable, as opposed to umbrellas glued externally to the fabric.

C. Grommet

The grommet links and secures both fabric layers together (D, E), forming a solid composite while providing a secure anchoring point for the umbrella.

D. U Band

The U Band prevents distal pistoning while allowing for the expansion of bulbous residual limbs.

E. External Fabric

The external fabric layer protects all internal elements and offers limited vertical stretch to provide a greater level of control and stability of the prosthesis.

F. Encapsulation

The aluminum encapsulation provides a highly durable attachment point to the prosthetic.

Learn more about our distal ends.

Take a deep dive into each distal end!






HOW ANTIOXIDANTS WORK

ALPS® EasyLiner® contains an antioxidant that is an effective scavenger of free radicals and hydroxyl groups.

Beneficial effects to the skin of patients have been reported. These stories have been provided here as anectodal evidence of the benefits that some patients have reported.





DIABETIC AMPUTEE

Case history provided by: Jim McElhiney, CPO, Nashville (USA).

This patient is a severe diabetic, insulin dependent with compromised vascularization. A local hospital called Jim McElhiney, CPO in Nashville (USA), to evaluate this patient.

The first two pictures have been taken prior to the fitting by Jim McElhiney. The patient had been out of her prosthesis for 3 weeks and was being treated with wet and dry dressing. She was then fitted with ALPS® EasyLiner® (Cushion Liner) and ALPS® Easysleeve (Sealing Sleeve), creating a suction suspension system.





Before the fitting with ALPS® EasyLiner®.





One week post initial fitting with ALPS® EasyLiner®.





Three weeks post initial fitting with ALPS® EasyLiner®.



PATIENT WITH DERMATITIS

Case history provided by: Wayne Sadler, CPO, Kentucky (USA).

This patient has been diagnosed with Contact Dermatitis. He followed the simple cleaning procedures as reccommended by the manufacturer.

When that wasn't enough, he started using our product, these are the results.





A DIFFICULT RESIDUAL LIMB

Case history provided by: Palgi Orthopaedic, Israel.





Before, when wearing a silicone liner.






After, when fitted with ALPS® EasyLiner® during 4 weeks.







ProFit Liner

Patent No. 1,076,089

ALPS® Pro-Fit is a prosthetic liner made with OptiGel that ensures excellent conformability to the residual limb. At the same time, its new and improved distal end encapsulation allows the user to experience a highly secure fit, resulting in enhanced control, comfort, and safety.



MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS

















OptiGel

Gray Fabric

K2 V K3 V K4 V

6 mm

Key Benefits



Made with OptiGel that has a smoother and softer texture than other liners on the market, ensuring gentle skin contact and a snug fit around the residual limb.



Features a new low-profile aluminum umbrella, allowing an enhanced interaction between the liner and limb with added distal cup flexibility to maximize user comfort.



Available in two uniform thicknesses, 3mm or 6mm, both of which feature thickened gel distally to enhance cushioning and comfort for the user.



Sleek, gray outer fabric affords superior wear and tear resistance while limiting pistoning due to its low vertical stretch capability.




ProFit Liner PFGDT/PFGFR. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	Locking 3mm or 6mm uniform	Cushion 3mm or 6mm uniform
16-19 cm	PFGDT16-(3 or 6)	PFGFR16-(3 or 6)
20-23 cm	PFGDT20-(3 or 6)	PFGFR20-(3 or 6)
24-25cm	PFGDT24-(3 or 6)	PFGFR24-(3 or 6)
26-27 cm	PFGDT26-(3 or 6)	PFGFR26-(3 or 6)
28-31 cm	PFGDT28-(3 or 6)	PFGFR28-(3 or 6)
32-34 cm	PFGDT32-(3 or 6)	PFGFR32-(3 or 6)
35-37 cm	PFGDT35-(3 or 6)	PFGFR35-(3 or 6)
38-43 cm	PFGDT38-(3 or 6)	PFGFR38-(3 or 6)
44-53 cm	PFGDT44-(3 or 6)	PFGFR44-(3 or 6)

Cross section








PFGDT

Locking

U-SHAPED REINFORCEMENT

Wider distal ends that are proportionate to the size of the liner



ALPS® Liner Technology





SoftSil® Liner

Patent No. 1,076,089, 11,602,445

ALPS® SoftSil Liner delivers the same medicalgrade qualities as our popular Silicone Pro but with enhanced features. Made with OptiSil, our new, soft, white silicone, and gray knitted fabric, the liner effectively offers a high degree of control and excellent freedom of movement without sacrificing a comfortable suspension in the user's prostheses.



MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS















OptiSil

Gray Fabric

K2 V K3 V K4 V

2 mm

Key Benefits



Medical grade, white silicone with 1/5 the friction of regular silicone offers superior comfort



Increased reduction of vertical stretch provides users with greater control and freedom of movement



Reduced pistoning



Skin friendly



Thin, lightweight & cooling




SoftSil Liner SSGSL/SSGSC. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	Locking 2.0 mm uniform	Cushion 2.0 mm uniform
20-21 cm	SSGSL-20	SSGSC-20
22-23 cm	SSGSL-22	SSGSC-22
24-25 cm	SSGSL-24	SSGSC-24
26-27 cm	SSGSL-26	SSGSC-26
28-29 cm	SSGSL-28	SSGSC-28
30-31 cm	SSGSL-30	SSGSC-30
32-34 cm	SSGSL-32	SSGSC-32
35-37 cm	SSGSL-35	SSGSC-35
38-40 cm	SSGSL-38	SSGSC-38
41-43 cm	SSGSL-41	SSGSC-41
44-47 cm	SSGSL-44	SSGSC-44

Cross section








SSGSL

Locking

U-SHAPED REINFORCEMENT

Wider distal ends that are proportionate to the size of the liner



ALPS® Liner Technology





UltraSeal Silicone Liner

Patent No. 1,076,089

The ALPS® UltraSeal Liner features an integrated band that creates an airtight seal as the user steps into the socket. The liner's smooth inner wall, constructed with OptiSil, a new soft, white silicone, easily conforms to and cushions the residual limb. The silicone, combined with a strongly knitted fabric covering, creates excellent durability and a secure fit.



MATERIAL





THICKNESS















OptiSil

Gray Fabric

K1 🗸

K3 < K4 <

2 mm

Key Benefits



Provides a comfortable cushion between the socket and residual limb and optimizes suspension while offering ultimate user convenience for daily use.



Outer knitted fabric combines strength and durability while allowing differentiated stretch and comfortable elasticity.



One sealing ring to ensure a secure hold.

Available at three different heights: Ultra Low, Low, and High



Medical-grade white silicone that contains 1/5 the friction of regular silicone



Skin friendly, lightweight, and cooling



UltraSeal Silicone Liner USHSL/USHSC/USLSL/USLSC/USUSL/USUSC. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	High Ring Locking 2mm uniform	High Ring Cushion 2mm uniform	Low Ring Locking 2mm uniform	Low Ring Cushion 2mm uniform	Ultra Low Ring Locking 2mm uniform	Ultra Low Ring Cushion 2mm uniform
20-21 cm	USHSL20	USHSC20	USLSL20	USLSC20	USUSL20	USUSC20
22-23 cm	USHSL22	USHSC22	USLSL22	USLSC22	USUSL22	USUSC22
24-25 cm	USHSL24	USHSC24	USLSL24	USLSC24	USUSL24	USUSC24
26-27 cm	USHSL26	USHSC26	USLSL26	USLSC26	USUSL26	USUSC26
28-29 cm	USHSL28	USHSC28	USLSL28	USLSC28	USUSL28	USUSC28
30-31 cm	USHSL30	USHSC30	USLSL30	USLSC30	USUSL30	USUSC30
32-34 cm	USHSL32	USHSC32	USLSL32	USLSC32	USUSL32	USUSC32
35-37 cm	USHSL35	USHSC35	USLSL35	USLSC35	USUSL35	USUSC35
38-40 cm	USHSL38	USHSC38	USLSL38	USLSC38	USUSL38	USUSC38
41-43 cm	USHSL41	USHSC41	USLSL41	USLSC41	USUSL41	USUSC41
44-47 cm	USHSL44	USHSC44	USLSL44	USLSC44	USUSL44	USUSC44

USHSL/USHSC: Ring begins 7in from the center of the distal end.

USLSL/USLSC: Ring begins 5in from the center of the distal end.


USUSL/USUSC: Ring begins 3in from the center of the distal end.

Cross section



USHSC USLSC USUSC

Hypobaric






USHSL USLSL USUSL

Locking

U-SHAPED REINFORCEMENT

Wider distal ends that are proportionate to the size of the liner



ALPS® Liner Technology





UltraSeal Gel Liner

Patent No. 1.076.089

Introducing the revolutionary UltraSeal gel liner—a groundbreaking product that seamlessly integrates advanced seal-in technology with the superior benefits of gel.

With an integrated band available in three different heights, users will experience an airtight seal upon socket entry. When combined with OptiGel, this liner offers exceptional comfort and excellent conformity to the residual limb, ensuring a secure and snug fit that enhances control, comfort, and safety.



MATERIAL



ACTIVITY LEVEL

THICKNESS















OptiGel

Gray Fabric

K1 🗸 K3 < K4 <

3 mm

Key Benefits



OptiGel has a smooth, soft texture that ensures gentle skin contact and a tailored fit around the residual limb.



Outer knitted fabric combines strength and durability while allowing differentiated stretch and comfortable elasticity.



The locking feature ensures users are provided with a suspension that is secure, firm, and comfortable.



One sealing ring to ensure a secure hold.

Available at three different heights: Ultra Low, Low, and High



UltraSeal Gel Liner USHGL/USHGC/USLGL/USLGC/USUGL/USUGC. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	High Ring Locking 3mm uniform	High Ring Cushion 3mm uniform	Low Ring Locking 3mm uniform	Low Ring Cushion 3mm uniform	Ultra Low Ring Locking 3mm uniform	Ultra Low Ring Cushion 3mm uniform
20-23 cm	USHGL20-3	USHGC20-3	USLGL20-3	USLGC20-3	USUGC20-3	USUGC20-3
24-25 cm	USHGL24-3	USHGC24-3	USLGL24-3	USLGC24-3	USUGC24-3	USUGC24-3
26-27 cm	USHGL26-3	USHGC26-3	USLGL26-3	USLGC26-3	USUGC26-3	USUGC26-3
28-31 cm	USHGL28-3	USHGC28-3	USLGL28-3	USLGC28-3	USUGC28-3	USUGC28-3
32-34 cm	USHGL32-3	USHGC32-3	USLGL32-3	USLGC32-3	USUGC32-3	USUGC32-3
35-37 cm	USHGL35-3	USHGC35-3	USLGL35-3	USLGC35-3	USUGC35-3	USUGC35-3
38-43 cm	USHGL38-3	USHGC38-3	USLGL38-3	USLGC38-3	USUGC38-3	USUGC38-3
44-47 cm	USHGL44-3	USHGC44-3	USLGL44-3	USLGC44-3	USUGC44-3	USUGC44-3

USHGL/USHGC: Ring begins 7in from the center of the distal end.

USLGL/USLGC: Ring begins 5in from the center of the distal end.


USUGL/USUGC: Ring begins 3in from the center of the distal end.

Cross section



USHGC USLGC USUGC

Hypobaric






USHGL USLGL USUGL

Locking

U-SHAPED REINFORCEMENT

Wider distal ends that are proportionate to the size of the liner



ALPS® Liner Technology





EasyLiner®

Patent No. 6,454,812* | 7,709,570 | 7,709,570

ALPS® EasyLiner® contains a prolonged-released antioxidant, an effective scavenger of free radicals and hydroxyl groups, with beneficial effects on the skin.

It provides the softness, strength, and cushioning from its use of ALPS® EasyGel, relieving the residual limb from the stress created by vertical and shear forces. It delicately conforms to the shape of the residual limb without restricting blood flow. It can be easily donned even by geriatric patients and users with reduced manual dexterity.



MATERIAL **FABRIC ACTIVITY LEVEL THICKNESS**



EasyGel

















Beige Fabric

3 mm

6 mm 6/3 mm distal proximal

Key Benefits



Indicated for delicate skin and diabetic patients thanks to the antioxidants contained in the gel



Protects the residual limb thanks to the high capacity of the EasyGel to absorb vertical and shear forces



Further relieves the skin from stress generated by shear forces thanks to its low coefficient of static friction



Delicately conforms to the shape of the residual limb and does not restrict blood flow



Easy to don also for patients with reduced manual dexterity




EasyLiner® ELDT/ELFR. Sizing Chart.

Measure the circumference of the stump at 6 cm above the distal end.


Measured size (cm)	Locking 3mm, 6mm, 6/3mm or 9mm	Cushion 3mm, 6mm, 6/3mm or 9mm
10-15 cm	N.A.	ELFR10-(3 or 6)
16-19 cm	ELDT16-(3, 6, or 6/3)	ELFR16-(3, 6, or 6/3)
20-23 cm	ELDT20-(3, 6, or 6/3)	ELFR20-(3, 6, or 6/3)
24-25 cm	ELDT24-(3, 6, 6/3, or 9)	ELFR24-(3, 6, 6/3, or 9)
26-27 cm	ELDT26-(3, 6, or 6/3)	ELRF26-(3, 6, or 6/3)
28-31 cm	ELDT28-(3, 6, 6/3, or 9)	ELRF28-(3, 6, 6/3, or 9)
32-34 cm	ELDT32-(3, 6, or 6/3)	ELFR32-(3, 6, or 6/3)
35-37 cm	ELDT35-(3 or 6)	ELFR35-(3 or 6)
38-43 cm	ELDT38-(3, 6, or 6/3)	ELFR38-(3, 6, or 6/3)
44-53 cm	ELDT44-(3, 6, or 6/3)	ELFR44-(3, 6, or 6/3)

Cross section









ALPS® Liner Technology





Silicone PRO Liner

Patent No. 5,507,834* | 6,454,812*

The Silicone PRO liner is manufactured with medical grade hypoallergenic silicone. The surface of the silicone features a low coefficient of static friction (1/5 of regular silicone). Without any additives, it is a hypoallergenic liner, very delicate on the skin. The liner is extremely lightweight and promotes the skin's natural thermoregulation.

The locking version features an innovative U matrix to limit vertical stretch while not constraining the residual limb.



MATERIAL



Silicone

FABRIC



Seamless Knitted Fabric **ACTIVITY LEVEL**



K1 🗸









/

THICKNESS



2 mm

Key Benefits



Does not constrict the residual limb and minimizes the pistoning effect thanks to the innovative U matrix (locking version)



Great freedom of movement thanks to the 2mm posterior thickness



 $Indicated \ for \ delicate \ skin \ since \ it \ is \ made \ of \ medical \ grade \ hypoallergenic \ silicone$



Minimizes friction to the residual limb due to its low coefficient of friction



Helps the thermoregulation of the residual limb



Easy to clean



Silicone PRO Liner LSC/LSL. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.


Measured size (cm)	Locking 2.0 mm uniform	Cushion 2.0 mm uniform
20-22 cm	LSL3.0-20	LSC3.0-20
22-24 cm	LSL3.0-22	LSC3.0-22
24-26 cm	LSL3.0-24	LSC3.0-24
26-28 cm	LSL3.0-26	LSC3.0-26
28-30 cm	LSL3.0-28	LSC3.0-28
30-32 cm	LSL3.0-30	LSC3.0-30
32-35 cm	LSL3.0-32	LSC3.0-32
35-38 cm	LSL3.0-35	LSC3.0-35
38-41 cm	LSL3.0-38	LSC3.0-38
41-44 cm	LSL3.0-41	LSC3.0-41
44-47 cm	LSL3.0-44	LSC3.0-44

Cross section



LSC

Cushion






LSL

Locking

U-SHAPED REINFORCEMENT

Wider distal ends that are proportionate to the size of the liner



ALPS® Liner Technology





Superior Performance Liner

Patent No. 6,454,812*

ALPS® Superior Performance prosthetic liner combines a limited vertical stretch fabric with our Grip Gel.

Grip Gel, with its softness, strength, and cushioning provides the greatest possible protection to the residual limb. The fabric, with its limited vertical elongation, improves control and stability of the prosthetic leg by reducing pistoning. The result is a comfortable liner, indicated for patients with low and medium activity levels.




MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS

















GripGel

Black Fabric

K1 🗸

K2 ✓ K3 ✓

3 mm 6 mm

Key Benefits



Protects the residual limb thanks to the high capacity of the Grip Gel to absorb vertical and shear forces



Improves control of the prosthetic leg thanks to the limited vertical stretch of the fabric



Reduces pistoning thanks to the fabric's limited vertical stretch and the distal matrix (locking version)



Safe suspension of the prosthetic leg thanks to the high level of adhesion of the Grip Gel




Superior Performance Liner SPDT/SPFR. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.


Measured size (cm)	Locking 3mm or 6mm uniform	Cushion 3mm or 6mm uniform
16-19 cm	SPDT16-(3 or 6)	SPFR16-(3 or 6)
20-23 cm	SPDT20-(3 or 6)	SPFR20-(3 or 6)
24-25 cm	SPDT24-(3 or 6)	SPFR24-(3 or 6)
26-27 cm	SPDT26-(3 or 6)	SPFR26-(3 or 6)
28-31 cm	SPDT28-(3 or 6)	SPFR28-(3 or 6)
32-34 cm	SPDT32-(3 or 6)	SPFR32-(3 or 6)
35-37 cm	SPDT35-(3 or 6)	SPFR35-(3 or 6)
38-43 cm	SPDT38-(3 or 6)	SPFR38-(3 or 6)
44-53 cm	SPDT44-(3 or 6)	SPFR44-(3 or 6)

Cross section









ALPS® Liner Technology





Superior Performance HD Liner

Patent No. 6,454,812*

ALPS® Superior Performance prosthetic liner, in its HD Gel version, provides the patient with comfort and control.

The HD Gel provides relieves the residual limb from vertical and shear forces; at the same time, it makes the liner more compact. The fabric, with its limited vertical elongation, further limits the vertical pistoning.



MATERIAL



ACTIVITY LEVEL

THICKNESS

















High Density Gel

Black Fabric

K2 V K3 V K4 V

6 mm

Key Benefits



Reduces pistoning thanks to the fabric's limited vertical stretch and the distal matrix (locking version)



Minimizes the movement of soft tissue and provides a high level of control of the prosthetic leg by combining the properties of the High Density Gel with those of the limited vertical stretch fabric



Protects the residual limb thanks to the high capacity of the HD Gel to absorb vertical and shear forces




Superior Performance Liner SPDTHD/SPFRHD. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	Locking 3mm or 6mm uniform	Cushion 3mm or 6mm uniform
16-19 cm	SPDTHD16-(3 or 6)	SPFRHD16-(3 or 6)
20-23 cm	SPDTHD20-(3 or 6)	SPFRHD20-(3 or 6)
24-25 cm	SPDTHD24-(3 or 6)	SPFRHD24-(3 or 6)
26-27 cm	SPDTHD26-(3 or 6)	SPFRHD26-(3 or 6)
28-31 cm	SPDTHD28-(3 or 6)	SPFRHD28-(3 or 6)
32-34 cm	SPDTHD32-(3 or 6)	SPFRHD32-(3 or 6)
35-37 cm	SPDTHD35-(3 or 6)	SPFRHD35-(3 or 6)
38-43 cm	SPDTHD38-(3 or 6)	SPFRHD38-(3 or 6)
44-53 cm	SPDTHD44-(3 or 6)	SPFRHD44-(3 or 6)

Cross section









ALPS® Liner Technology





General Purpose Liner

Patent No. 6,454,812*

The General Purpose liner contains the softness, strength, and cushioning of our Grip Gel, relieving the residual limb from the stress generated by vertical and shear forces The superior elongation of this prosthetic liner means that blood flow is not restricted and it allows the liner to delicately conform to the shape of the residual limb. It is easily donned even by geriatric patients or those with reduced manual dexterity.



MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS

















GripGel

Beige Fabric

K1 ✓

K2 V K3 V K4 V

6 mm

Key Benefits



Delicately conforms to the shape of the residual limb and does not restrict blood flow



Easy to don also for patients with reduced manual dexterity



Protects the residual limb thanks to the high capacity of the Grip Gel to absorb vertical and shear forces



Safe suspension of the prosthesis due to the superior adhesive properties of the Grip Gel



General Purpose Liner GPDT/GPFR. Sizing Chart.


Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	Locking 3mm or 6mm uniform	Cushion 3mm or 6mm uniform
16-19 cm	GPDT16-(3 or 6)	GPFR16-(3 or 6)
20-23 cm	GPDT20-(3 or 6)	GPFR20-(3 or 6)
24-25 cm	GPDT24-(3 or 6)	GPFR24-(3 or 6)
26-27 cm	GPDT26-(3 or 6)	GPRF26-(3 or 6)
28-31 cm	GPDT28-(3 or 6)	GPRF28-(3 or 6)
32-34 cm	GPDT32-(3 or 6)	GPFR32-(3 or 6)
35-37 cm	GPDT35-(3 or 6)	GPFR35-(3 or 6)
38-43 cm	GPDT38-(3 or 6)	GPFR38-(3 or 6)
44-53 cm	GPDT44-(3 or 6)	GPFR44-(3 or 6)
44-70 cm*	GPDT-XL*	N.A.


^{*}GPDT-XL comes in 3mm thickness only. It features a larger umbrella and a U-Shaped matrix. Ideal for bariatric patients.

Cross section









ALPS® Liner Technology





Eco™ Liner

Patent No. 6,454,812* | 7,709,570 | 10,966,845 | 7,709,570

ALPS® Eco™ Liner eatures an innovative knitted fabric with areas of varying elasticity to improve the control of the prosthetic leg. The HD Gel helps keep the residual limb more compact and provides enhanced stability to the prosthetic leg, at the same time ensuring a high level of protection from vertical and shear forces. The locking version features an innovative U matrix to limit the vertical stretch without constraining the residual limb.



MATERIAL



High Density Gel

FABRIC



Seamless Knitted Fabric

ACTIVITY LEVEL



K1 🗸









THICKNESS





6 mm

Key Benefits



Protects the residual limb thanks to the high capacity of the HD Gel to absorb vertical and shear forces



Ensures maximum control of the prosthetic leg and minimizes soft tissue movement thanks to the High Density Gel



Improved comfort due to the effect of the innovative seamless knitted fabric providing areas of varying levels of elasticity



Does not constrict the residual limb and minimizes the pistoning effect thanks to the innovative U matrix (locking version)




Eco[™] Liner ECDTHD/ECFRHD. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	Locking 3mm or 6mm uniform	Cushion 3mm or 6mm uniform
20-23 cm	ECDTHD20-(3 or 6)	ECFRHD20-(3 or 6)
24-25 cm	ECDTHD24-(3 or 6)	ECFRHD24-(3 or 6)
26-27 cm	ECDTHD26-(3 or 6)	ECFRHD26-(3 or 6)
28-31 cm	ECDTHD28-(3 or 6)	ECFRHD28-(3 or 6)
32-34 cm	ECDTHD32-(3 or 6)	ECFRHD32-(3 or 6)
35-37 cm	ECDTHD35-(3 or 6)	ECFRHD35-(3 or 6)
38-43 cm	ECDTHD38-(3 or 6)	ECFRHD38-(3 or 6)
44-53 cm	ECDTHD44-(3 or 6)	ECFRHD44-(3 or 6)

Cross section








ECDTHD

Locking

U-SHAPED REINFORCEMENT

Wider distal ends that are proportionate to the size of the liner



ALPS® Liner Technology





Smart Seal™ Liner

Patent No. 6,454,812* | 8,852,291 | 9,216,099 | 11,179,252 | 8,852,291 9,216,099 | 11,179,252

ALPS® Smart Seal liner features integrated, silicone coated sealing bands that create a safe interface.

The bands stretch against the socket wall, preventing the liner from slipping.

The inner HD Gel layer, thanks to its low modulus of elasticity, easily conforms to the residual limb without restricting the blood flow.

It is available in two configurations with different sealing band heights.



MATERIAL



High Density Gel

FABRIC



Seamless Knitted Fabric

ACTIVITY LEVEL









THICKNESS





3 mm 6 mm

Key Benefits



Ensures a high level of control of the prosthetic leg and minimizes soft tissue movement thanks to the High Density Gel



Non-constrictive hypobaric suspension thanks to its integrated, silicone coated sealing bands



Indicated both for transtibial and transfemoral users since the sealing bands are available in two different configurations



Maximum stability of the liner within the socket due to the High Density Gel and the use of a limited vertical stretch fabric in the distal area



Smart Seal[™] Liner SAFR/SADT/SBFR. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	SAFR - Cushion 3mm or 6mm uniform
26-27 cm	SAFR26-(3 or 6)
28-31 cm	SAFR28-(3 or 6)
32-34 cm	SAFR32-(3 or 6)
35-37 cm	SAFR35-(3 or 6)
38-43 cm	SAFR38-(3 or 6)

Measured size (cm)	SBFR - Cushion 3mm or 6mm uniform
20-23 cm	SBFR20-(3 or 6)
24-25 cm	SBFR24-(3 or 6)
26-27 cm	SBFR26-(3 or 6)
28-31 cm	SBFR28-(3 or 6)

SAFR

Rings begin 12.5 cm from the center of the distal end.

SBFR

Rings begin 7.5 cm from the center of the distal end.


SADT

Sizes 20-28: rings begin 7.5 cm from the center of the distal end. Sizes 32-38: rings begin 12.7 cm from the center of the distal end.

Measured size (cm)	SADT - Locking 3mm or 6mm uniform
20-23 cm	SADT20-(3 or 6)
24-25 cm	SADT24-(3 or 6)
26-27 cm	SADT26-(3 or 6)
28-31 cm	SADT28-(3 or 6)
32-37 cm	SADT32-(3 or 6)
38-43 cm	SADT38-(3 or 6)

Cross section








SADT

Locking

U-SHAPED REINFORCEMENT

Wider distal ends that are proportionate to the size of the liner



ALPS® Liner Technology





ExtremeTM Liner

Patent No. 6,454,812*

The ALPS® Extreme prosthetic liner combines softness, strength, and cushioning from our Grip Gel with a fabric with limited vertical stretch. It improves the control of the prosthetic leg and reduces pistoning, while still providing superior protection of the residual limb. The superior adhesion of the Grip Gel allows the safe suspension of the prosthetic leg. The result is a comfortable and versatile liner, indicated for patients with low and medium activity levels.




MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS

















6 mm

GripGel

Brown Fabric

K1 ✓

K2 ✓ K3 ✓

3 mm

Key Benefits



Protects the residual limb thanks to the high capacity of the Grip Gel to absorb vertical and shear forces



Improves control of the prosthetic leg thanks to the limited vertical stretch of the fabric



Reduces pistoning thanks to the fabric's limited vertical stretch and the distal matrix (locking version)



Safe suspension of the prosthesis due to the superior adhesive properties of the Grip Gel




Extreme™ Liner AKDT/AKFR. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.


Measured size (cm)	Locking 3mm or 6mm uniform	Cushion 3mm or 6mm uniform
16-19 cm	AKDT16-(3 or 6)	AKFR16-(3 or 6)
20-23 cm	AKDT20-(3 or 6)	AKFR20-(3 or 6)
24-25 cm	AKDT24-(3 or 6)	AKFR24-(3 or 6)
26-27 cm	AKDT26-(3 or 6)	AKFR26-(3 or 6)
28-31 cm	AKDT28-(3 or 6)	AKFR28-(3 or 6)
32-34 cm	AKDT32-(3 or 6)	AKFR32-(3 or 6)
35-37 cm	AKDT35-(3 or 6)	AKFR35-(3 or 6)
38-43 cm	AKDT38-(3 or 6)	AKFR38-(3 or 6)
44-53 cm	AKDT44-(3 or 6)	AKFR44-(3 or 6)

Cross section









ALPS® Liner Technology





Extreme™ HD Liner

Patent No. 6,454,812*

ALPS® Extreme™ prosthetic liner, in its HD Gel version, provides the patient with comfort and control. The HD Gel relieves the residual limb from vertical and shear forces; at the same time, it makes the liner more compact. The fabric, with its limited vertical elongation, further limits vertical pistoning. This makes it an ideal liner for the most active users who are looking to maintain maximum control of the prosthetic leg without compromising comfort.



MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS

















High Density Gel

Brown Fabric

K2 V K3 V K4 V

6 mm

Key Benefits



Reduces pistoning thanks to the fabric's limited vertical stretch and the distal matrix (locking version)



Minimizes the movement of soft tissue and provides a high level of control of the prosthetic leg by combining the properties of the High Density Gel with those of the limited vertical stretch fabric



Protects the residual limb thanks to the high capacity of the HD Gel to absorb vertical and shear forces




Extreme™ Liner AKDTHD/AKFRHD. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.


Measured size (cm)	Locking 3mm or 6mm uniform	Cushion 3mm or 6mm uniform
16-19 cm	AKDTHD16-(3 or 6)	AKFRHD16-(3 or 6)
20-23 cm	AKDTHD20-(3 or 6)	AKFRHD20-(3 or 6)
24-25 cm	AKDTHD24-(3 or 6)	AKFRHD24-(3 or 6)
26-27 cm	AKDTHD26-(3 or 6)	AKFRHD26-(3 or 6)
28-31 cm	AKDTHD28-(3 or 6)	AKFRHD28-(3 or 6)
32-34 cm	AKDTHD32-(3 or 6)	AKFRHD32-(3 or 6)
35-37 cm	AKDTHD35-(3 or 6)	AKFRHD35-(3 or 6)
38-43 cm	AKDTHD38-(3 or 6)	AKFRHD38-(3 or 6)
44-53 cm	AKDTHD44-(3 or 6)	AKFRHD44-(3 or 6)

Cross section









ALPS® Liner Technology





EasyFlex[™] Liner

Patent No. 6,454,812* | 7,709,570 | 8,226,732 | 8,246,694 | 8,349,021

ALPS® EasyFlex™ prosthetic liner allows users ease of knee flexion thanks to an unlimited restriction of stretch anteriorly. Posteriorly the fabric, with its limited vertical stretch, improves the control of the prosthetic leg. The prolonged-released antioxidant contained in the EasyGel is an effective scavenger of free radicals and hydroxyl groups, providing beneficial effects on the skin. The EasyFlex[™] Liner is designed for use with a pin and lock suspension system.



MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS





















EasyGel

Anterior Posterior Beige Brown Fabric Fabric

K1 🗸

K3 🗸 K4 🗸 K2 🗸

3 mm 6 mm

Key Benefits



Indicated for delicate skin and diabetic patients thanks to the antioxidants contained in the gel



Protects the residual limb thanks to the high capacity of the EasyGel to absorb vertical and shear forces



Further relieves the skin from stress generated by shear forces thanks to its low coefficient of static friction



Allows an easy knee flexion and improves the control of the prosthetic leg, since it has been engineered using two types of fabric with different elasticity




EasyFlexTM Liner EZF. Sizing Chart. Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	Locking 3mm or 6mm uniform
16-19 cm	EZF16-(3 or 6)
20-23 cm	EZF20-(3 or 6)
24-25 cm	EZF24-(3 or 6)
26-27 cm	EZF26-(3 or 6)
28-31 cm	EZF28-(3 or 6)
32-34 cm	EZF32-(3 or 6)
35-37 cm	EZF35-(3 or 6)
38-43 cm	EZF38-(3 or 6)
44-53 cm	EZF44-(3 or 6)

Cross section





ALPS® Liner Technology





EasyLiner® SuperStretch

Patent No.7,709,570

ALPS® EasyLiner® SuperStretch is a soft, robust and versatile prosthetic liner made of EasyGel without the restriction of a fabric cover.

It allows complete freedom of movement and conforms to irregularly shaped residual limbs without restricting blood flow. Due to the low coefficient of friction of the EasyGel, it relieves the skin from the stress created by shear forces. The prolonged-released antioxidant has beneficial effects on the skin.



MATERIAL

ACTIVITY LEVEL

THICKNESS



EasyGel





K1✓ K2✓ K3









6 mm

Key Benefits



Indicated for delicate skin and diabetic patients thanks to the antioxidants contained in the gel



Protects the residual limb thanks to the high capacity of the EasyGel to absorb vertical and shear forces



Further relieves the skin from stress generated by shear forces thanks to its low coefficient of static friction



Made of EasyGel with no fabric cover, it allows a superior fit to irregular shapes without constricting the residual limb




EasyLiner® SuperStretch ELPX. Sizing Chart. Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	Cushion 6mm uniform
10-15 cm	ELPX-10
16-19 cm	ELPX-16
20-23 cm	ELPX-20
24-25 cm	ELPX-24
26-27 cm	ELPX-26
28-31 cm	ELPX-28
32-37 cm	ELPX-32
38-43 cm	ELPX-38

Cross section





ALPS® Liner Technology





Anterior Posterior Tapered Liner

Patent No.6,454,812* | 8,226,732 | 8,246,694 | 8,349,021 | 8,226,732

ALPS® AP Tapered Liner is gradually tapered from the thicker anterior to the thinner posterior side to provide freedom of movement to transtibial users.

The 6 mm frontal wall ensures a high level of protection to bony areas. By combining a high density gel with a limited vertical stretch fabric, this liner minimizes soft tissue movement and reduces vertical pistoning, offering a high level of control of the prosthetic leg.



MATERIAL



ACTIVITY LEVEL

THICKNESS

















High Density Gel

Black Fabric

K2 V K3 V K4 V

6/3 mm Anterior Posterior Taper

Key Benefits



Reduces pistoning thanks to the fabric's limited vertical stretch and the distal matrix (locking version)



Protects the residual limb thanks to the high capacity of the HD Gel to absorb vertical and shear forces



Protects bony areas with its 6 mm anterior wall thickness



Great freedom of movement thanks to the 3 mm posterior thickness



Minimizes soft tissue movement and provides superior control of the prosthetic leg by combining High Density Gel with a limited vertical stretch fabric




Anterior Posterior Tapered Liner APDT/APFR. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.


Measured size (cm)	Locking 6/3mm tapered	Cushion 6/3mm tapered
20-23 cm	APDT20-TP	APFR20-TP
24-25 cm	APDT24-TP	APFR24-TP
26-27 cm	APDT26-TP	APFR26-TP
28-31 cm	APDT28-TP	APFR28-TP
32-37 cm	APDT32-TP	APFR32-TP

Cross section









ALPS® Liner Technology





PSI Liner

Patent No.7,709,570 | 8,357,206 | 8,394,150

The PSI Liner improves suspension by sealing onto the sleeves of a transtibial prosthesis. The EasyGel provides comfort and durability, relieving the skin from the stress created by vertical and shear forces. The antioxidant contained in the gel has beneficial effects on the skin.

The PSI Liner can be used in conjunction with current valve suction or vacuum systems being used on the market today.



MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS









K2 🗸









EasyGel

Black Fabric

1~

K

K/i

6 mm 6/3 mm
distal proximal

Key Benefits



Indicated for delicate skin and diabetic patients thanks to the antioxidants contained in the gel



Protects the residual limb thanks to the high capacity of the EasyGel to absorb vertical and shear forces



Further relieves the skin from stress generated by shear forces thanks to its low coefficient of static friction



The proximal area with no fabric cover seals onto the sleeve to create a secure suspension




Proximal Seal Interface Liner PSI. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	Cushion 6mm or 3/6mm tapered
16-19 cm	PSI16-(6 or 6/3)
20-23 cm	PSI20-(6 or 6/3)
24-25 cm	PSI24-(6 or 6/3)
26-27 cm	PSI26-(6 or 6/3)
28-31 cm	PSI28-(6 or 6/3)
32-37 cm	PSI32-(6 or 6/3)
38-43 cm	PSI38-(6 or 6/3)
44-53 cm	PS144-(6 or 6/3)

Cross section





ALPS® Liner Technology

 ${f 5}$ Elastomers. ${f 5}$ Fabrics. ${f 5}$ Suspensions.





Liberty™ Liner

Patent No. 6,544,292* | 6,454,812* | 6,361,568* | 8,382,852 | 8,808,394

The Liberty™ Liner allows the expulsion of air and humidity through a built-in valve integrated either on the pin or the reducer, promoting a healthy environment.

Used in conjunction with a nylon sheath, it presents a superior capacity for collecting and expelling humidity.

The HD gel improves the stability of the prosthetic leg and provides greater protection against vertical and shear forces. It is recommended to minimize the problem of excessive perspiration.



MATERIAL



High Density Gel

FABRIC



Seamless Knitted Fabric

ACTIVITY LEVEL



K1 🗸





K2 V K3 V K4 V







THICKNESS



6 mm

Key Benefits



Protects the residual limb thanks to the high capacity of the HD Gel to absorb vertical and shear forces



Promotes a healthy environment due to the distal valve that expels air and humidity between the residual limb and the liner



Ensures a high level of control of the prosthetic leg and minimizes soft tissue movement thanks to the High Density Gel



Since it is used in conjunction with a nylon sheath, it has a superior capacity for collecting and expelling humidity



Liberty[™] Liner II VSDT. Sizing Chart.

Measure the circumference of the residual limb at 6 cm above the distal end.

Measured size (cm)	Locking 3mm uniform	Locking 6mm uniform
20-23 cm	VSDT20-3	VSDT20-6
24-25 cm	VSDT24-3	VSDT24-6
26-27 cm	VSDT26-3	VSDT26-6
28-31 cm	VSDT28-3	VSDT28-6
32-37 cm	VSDT32-3	VSDT32-6
38-43 cm	VSDT38-3	VSDT38-6

Cross section




VSDT

Locking

U-SHAPED REINFORCEMENT

Wider distal ends that are proportionate to the size of the liner



ALPS® Liner Technology









Reinforced Flex Sleeve

The ALPS® Reinforced Flex prosthetic suspension sleeve features a 30° pre-flexion design and includes an additional layer of internal fabric for enhanced durability around the trim line of the user's socket. Additionally, the combination of OptiGel and gray knitted fabric provides maximum comfort and high performance, making the RFX® an excellent choice for active



MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS

















3 mm

OptiGel

Gray or Black

K2 V K3 V K4 V

Key Benefits



Extra layer of internal fabric provides a more durable option for high-activity users



Gray knitted fabric combined with OptiGel delivers maximum comfort, control, and high performance



30° pre-flexion to ease bending and reduce bunching behind the knee



Provides a secure suspension



Available in two color options (black and gray)



Available in one length that accommodates most users




Reinforced Flex Sleeve RFX®. Sizing Chart.

Measure the circumference of the leg at 10 cm above the center of the knee.

Measured size (cm)	3mm gray	3mm black
34-37 cm	RFX20-3G	RFX20-3BL
36-39 cm	RFX23-3G	RFX23-3BL
38-42 cm	RFX26-3G	RFX26-3BL
41-46 cm	RFX29-3G	RFX29-3BL
45-50 cm	RFX32-3G	RFX32-3BL

Cross section



ALPS® Liner Technology

 $\bf 5$ Elastomers. $\bf 5$ Fabrics. $\bf 5$ Suspensions.





Reinforced Sleeve

ALPS® Reinforced Sleeve is constructed with gray knitted fabric, OptiGel, and an extra layer of internal fabric to deliver superior durability. This combination works together to provide the maximum comfort, control, and high-performance capabilities required for higher activity-level users.



MATERIAL

FABRIC

ACTIVITY LEVEL

THICKNESS















OptiGel

Gray Fabric

K2 V K3 V K4 V

3 mm

Key Benefits



Extra layer of internal fabric provides increased durability



Gray knitted fabric's differential stretch attributes allow for maximum comfort and stability



Provides a secure suspension



OptiGel offers a great degree of control



The gel sticks to the patients skin without causing shear force and helps prevent the sleeve from rolling down the amputee's limb.




Reinforced Sleeve RFS. Sizing Chart.

Measure the circumference of the leg at 10 cm above the center of the knee.

Measured size (cm)	3mm gray
25-30 cm	RFS20-3G
31-34 cm	RFS24-3G
35-37 cm	RFS26-3G
38-40 cm	RFS28-3G
41-45 cm	RFS32-3G
46-60 cm	RFS38-3G
61-70 cm	RFS44-3G

Cross section



ALPS® Liner Technology

5 Elastomers. 5 Fabrics. 5 Suspensions.





Flex Sleeve

ALPS® Flex Sleeve (SFX) prosthetic suspension sleeve combines a seamless knitted fabric with a High Density Gel to provide maximum comfort and high performance.

The SFX provides a 30° pre-flexion which eases bending and reduces bunching behind the knee. It is indicated for medium and high activity users.



MATERIAL



High Density Gel

FABRIC



Seamless Knitted Fabric **ACTIVITY LEVEL**











THICKNESS



3 mm

Key Benefits



Minimizes compression and irritation proximally



Ease of knee flexion and reduced bunching in the popliteal



The High Density Gel ensures superior control even for more active users



Available in two color options (beige and gray)



Provides a secure suspension of the prosthetic leg



Recommended for waterproof prostheses




Flex Sleeve SFX. Sizing Chart.

Measure the circumference of the leg at 10 cm above the center of the knee.

Measured size (cm)	3mm beige	3mm gray
34-37 cm	SFX20-3B	SFX20-3G
36-39 cm	SFX23-3B	SFX23-3G
38-42 cm	SFX26-3B	SFX26-3G
41-46 cm	SFX29-3B	SFX29-3G
45-50 cm	SFX32-3B	SFX32-3G

Cross section



ALPS® Liner Technology

 $\bf 5$ Elastomers. $\bf 5$ Fabrics. $\bf 5$ Suspensions.





EasySleeve

Patent No. 7.709.570

ALPS® EasySleeve Reinforced Fabric prosthetic gel sleeve provides superior comfort thanks to its use of ALPS® EasyGel. The elasticity of the fabric allows the user a great deal of freedom of knee flexion. The EasySleeve also achieves secure suspension by sealing to the skin without restricting circulation. The lower coefficient of friction found with the EasyGel helps relieve the skin of shear forces.



MATERIAL



ACTIVITY LEVEL

THICKNESS







Beige Fabric













EasyGel

K2 V K3 V K4 V

6 mm

Key Benefits



Easy to don also for patients with reduced manual dexterity



Indicated for delicate and fragile skins thanks to the antioxidants contained in the gel



Minimizes compression and irritation proximally



Allows a high degree of freedom of knee flexion




Relieves the skin from stress generated by shear forces thanks to EasyGel's low coefficient of static friction



EasySleeve SFR. Sizing Chart. Measure the circumference of the leg at 10 cm above the center of the knee.

Measured size (cm)	3mm or 6mm uniform
16-19 cm	SFR10-(3 or 6)
20-24 cm	SFR16-(3 or 6)
25-30 cm	SFR20-(3 or 6)
31–34 cm	SFR24-(3 or 6)
35-37 cm	SFR26-(3 or 6)
38-40 cm	SFR28-(3 or 6)
41–45 cm	SFR32-(3 or 6)
46-60 cm	SFR38-(3 or 6)
61-70 cm	SFR44-(3 or 6)

Cross section



ALPS® Liner Technology

5 Elastomers. 5 Fabrics. 5 Suspensions.





EasySleeve SuperStretch

Patent No. 7,709,570

ALPS® EasySleeve SuperStretch is a prosthetic gel sleeve that allows freedom of knee flexion as well as providing secure suspension by sealing to the skin without restricting circulation.

By using a separate fabric cover, the EasyGel is able to move with the user, therefore reducing the effect of physical forces against the patella.



MATERIAL

ACTIVITY LEVEL

THICKNESS

















6 mm

Key Benefits



Indicated for delicate and fragile skins thanks to the antioxidants contained in the gel



Minimizes compression and irritation proximally



Allows a high degree of freedom of knee flexion



Relieves the skin from stress generated by shear forces thanks to EasyGel's low coefficient of static friction



A fabric cover (sold separately) protects the sleeve



EasySleeve SuperStretch SLE. Sizing Chart.

Measure the circumference of the leg at 10 cm above the center of the knee.

Measured size (cm)	6mm uniform	Cover size
16-19 cm	SLE-10	N.A.
20-24 cm	SLE-16	CSLE-16
25-30 cm	SLE-20	CSLE-20
31-34 cm	SLE-24	CSLE-24
35-37 cm	SLE-26	CSLE-26
38-40 cm	SLE-28	CSLE-28
41-45 cm	SLE-32	CSLE-32

Cross section



ALPS® Liner Technology

5 Elastomers. 5 Fabrics. 5 Suspensions.





Superior Suspension Sleeve

ALPS® Superior Suspension Sleeve provides superior comfort thanks to its use of ALPS® Gel.

The Gel adheres to the amputee's skin to achieve secure suspension without restricting circulation and without causing shear forces.

With its thinner profile, the SFB, is the ideal choice for patients who want a lightweight product with reduced bulk.



MATERIAL





High Density

GripGel

ACTIVITY LEVEL









THICKNESS



2 mm

Key Benefits



Minimizes compression and irritation proximally



The High Density Gel ensures superior control even for more active users



Provides a secure suspension of the prosthetic leg



Reduced bulk thanks to its 2 mm profile



The fabric's limited vertical stretch improves control of the prosthetic leg and allows for ease of knee flexion




Superior Suspension Sleeve SFB/SFBHD. Sizing Chart.

Measure the circumference of the leg at 10 cm above the center of the knee.

Measured size (cm)	2mm uniform
25-30 cm	SFB(HD)20-2
31-34 cm	SFB(HD)24-2
35-37 cm	SFB(HD)26-2
38-40 cm	SFB(HD)28-2
41-45 cm	SFB(HD)32-2

Cross section





5 Elastomers. **5** Fabrics. **5** Suspensions.





Superior Performance Sleeve

ALPS® Superior Performance Sleeve is a prosthetic gel sleeve constructed using a limited vertical stretch fabric that maintains a high horizontal elongation to provide the transtibial user with easy knee flexion and elevated control.

The Gel adheres to the patient's skin without creating shear force and helps prevent the sleeve from rolling down the residual limb.



MATERIAL







GripGel

ACTIVITY LEVEL











THICKNESS





6 mm

Key Benefits



Minimizes compression and irritation proximally



The High Density Gel ensures superior control even for more active users



Provides a secure suspension of the prosthetic leg



The fabric's limited vertical stretch improves control of the prosthetic leg and allows for ease of knee flexion




Superior Performance Sleeve SFP/SFPHD. Sizing Chart.

Measure the circumference of the leg at 10 cm above the center of the knee.

3mm or 6mm uniform
SFP(HD)16-(3 or 6)
SFP(HD)20-(3 or 6)
SFP(HD)24-(3 or 6)
SFP(HD)26-(3 or 6)
SFP(HD)28-(3 or 6)
SFP(HD)32-(3 or 6)
SFP(HD)38-(3 or 6)
SFP(HD)44-(3 or 6)

Cross section











Extreme Sleeve

ALPS® Extreme Sleeve is covered with a limited vertical stretch fabric that maintains a high horizontal elongation to provide the user with an easy knee flexion and elevated control.

The Grip Gel adheres to the patient's skin without creating shear force and helps prevent the sleeve from rolling down the residual limb.



MATERIAL

ACTIVITY LEVEL

THICKNESS















GripGel

K1 V K2 V K3 V K4 V

3 mm 6 mm

Key Benefits



Minimizes compression and irritation proximally



Provides a secure suspension of the prosthetic leg



The fabric's limited vertical stretch improves control of the prosthetic leg and allows for ease of knee flexion




Extreme Sleeve SFK. Sizing Chart.

Measure the circumference of the leg at 10 cm above the center of the knee.

3mm or 6mm uniform
SFK16-(3 or 6)
SFK20-(3 or 6)
SFK24-(3 or 6)
SFK26-(3 or 6)
SFK28-(3 or 6)
SFK32-(3 or 6)
SFK38-(3 or 6)
SFK44-(3 or 6)

Cross section





 ${f 5}$ Elastomers. ${f 5}$ Fabrics. ${f 5}$ Suspensions.





V.I.V.A. Sleeve

Patent No. 6,361,568*

ALPS® V.I.V.A. Sleeve is a prosthetic gel sleeve with a built-in valve to evacuate air from the socket to ensure superior security.

It provides suction suspension without the need to create a valve in the socket.

Once the sleeve is donned, the valve should be located proximal to the patella to create back pressure with knee flexion.



MATERIAL

ACTIVITY LEVEL

THICKNESS















GripGel

K1✓ K2✓ K3

6 mm

Key Benefits



Minimizes compression and irritation proximally



Allows a high degree of freedom of knee flexion



Provides a secure suspension of the prosthetic leg




Improves pneumatic suspension thanks to the built-in valve for evacuating air



V.I.V.A. Sleeve VIVA. Sizing Chart. Measure the circumference of the leg at 10 cm above the center of the knee.

Measured size (cm)	3mm or 6mm uniform
20-24 cm	VIVA16-(3 or 6)
25-30 cm	VIVA20-(3 or 6)
31-34 cm	VIVA24-(3 or 6)
35-37 cm	VIVA26-(3 or 6)
38-40 cm	VIVA28-(3 or 6)
41-45 cm	VIVA32-(3 or 6)

Cross section





5 Elastomers. 5 Fabrics. 5 Suspensions.





Anterior Posterior Tapered Sleeve

The ALPS Anterior Posterior Tapered Sleeve is gradually tapered from the anterior to the posterior to provide superior comfort. It features our black high performance fabric with HD gel to provide amputees with a great degree of control.



MATERIAL





ACTIVITY LEVEL





THICKNESS



High Density Gel

K1 × K2 × K3 × K4 ×

6/3mm AP Tapered

Key Benefits



Black fabric offers horizontal stretch with highly reduced vertical stretch



Allows a high degree of freedom of knee flexion



Provides a secure suspension of the prosthetic leg




HD Gel provides active amputees a great degree of control



Anterior Posterior Tapered Sleeve. APS. Sizing Chart. Measure the circumference of the leg at 10 cm above the center of the knee.

Measured size (cm)	6/3 mm tapered
20-23 cm	APS20-TP
24-25 cm	APS24-TP
26-27 cm	APS26-TP
28-31 cm	APS28-TP
32-37 cm	APS32-TP

Cross section





5 Elastomers. 5 Fabrics. 5 Suspensions.





EasySleeve Fabric Reinforced

ALPS® SFS has superior comfort provided by GripGel. The elasticity of the fabric allows the user excellent freedom of knee flexion. The SFS also achieves secure suspension by sealing with the skin without restricting circulation, while the GripGel sticks to the patients skin without causing shear force.



MATERIAL

ACTIVITY LEVEL

THICKNESS









K1 ✓ K2 ✓ K3











3 mm 6 mm

Key Benefits



Beige, high-performance fabric improves the sleeve's durability



Allows a high degree of freedom of knee flexion



Provides a secure suspension of the prosthetic leg



Minimizes compression and irritation proximally




The GripGel helps prevent the sleeve from rolling down the amputee's limb



EasySleeve Fabric Reinforced. SFS. Sizing Chart. Measure the circumference of the leg at 10 cm above the center of the knee.

Measured size (cm)	3mm or 6mm uniform
20-24 cm	SFS16-(3 or 6)
25-30 cm	SFS20-(3 or 6)
31-34 cm	SFS24-(3 or 6)
35-37 cm	SFS26-(3 or 6)
38-40 cm	SFS28-(3 or 6)
41-45 cm	SFS32-(3 or 6)
46-60 cm	SFS38-(3 or 6)
61-70 cm	SFS44-(3 or 6)

Cross section











Vacuum Integrated Pump

Patent No. 8,197,555 | 11,224,527 | 8,197,555 | 11,224,527

The ALPS® VIP® is an air expulsion integrated pump for elevated vacuum. Its use provides patients continuous total contact when donning the prosthesis. This is essential for an even weight distribution and to distribute the suspension of the prosthesis over the whole surface of the residual limb. Elevated vacuum created by the VIP® stabilizes the residual limb, overcoming all issues connected with the residual limb changing volume during the day.



Key Benefits



Efficiently deals with residual limb volume fluctuation which in turn leads to improved prosthesis control and comfort



Reduces the intensity of peaks in pressure on the patient's limbs, improving comfort and providing greater freedom of movement



Helps preserve residual limb skin health by improving circulation and consequently increasing limb oxygenation



Quiet, light, and does not affect the cosmetic appearance of the prosthesis



Easy to install and battery-free




Vacuum Integrated Pump VIP® Part Number and Optional Complementary Products.

Available in 2 different vacuum levels: HIGH VACUUM (Red) 10-13 in of mercury and LOW VACUUM (Green) 5-8 in of mercury. Choose vacuum level on the basis of activity level and skin condition.

Part name	Part #
VIP® LOW VACUUM (GREEN)	VIP - DL
VIP® HIGH VACUUM (RED)	VIP - DH
DUMMY	VE319 - D
SHORT DUMMY	VE319 - SD
MOUNTING PLATE	VE819 - P
LOCK (SMALL)	VIP-L1120-S
MINI VALVE	MV-009
MINI VALVE KIT	LMV-009
PIN (4.8cm)	VIP-92L
PIN (3.0cm)	VIP-92S

Cross section

With Mounting Plate



Without Mounting Plate





VIP® Lock

With ALPS® VIP® Lock, clinicians can quickly switch from pin & lock to vacuum suspension, or vice versa, with minor adjustments to the socket.



Key Benefits



Retrofits with ALPS® VIP®



Lightweight



Removable release pin




Perfect fit with ALPS® encapsulations



VIP® Lock Part Number and Specification

PART #	VIP-L1120-S		
DIMENSIONS	32x67 mm		
WEIGHT	74 g		
WEIGHT LIMIT	270 kg		
PIN (4.8cm)	VIP-92L		
PIN (3.0cm)	VIP-92S		





Lock with Pyramid

Patent No. 5,507,837*

The ALPS® Lock with Pyramid is an all-in-one lock body design which combines the lock and pyramid. The glass reinforced nylon lock body exhibits better wear characteristics than delrin or aluminum. The convenient one-way gear allows the user to push in or wind out the limb into the prosthetic socket. A unique feature of this lock is that if foreign material or clothing get caught in the gear, the housing can be removed, freeing the pin.



Key Benefits



For laminated sockets



All-in-one lock body design combining the lock and pyramid



Infinite pin adjustment reduces play



The housing can be removed, freeing the pin




100 kg (220 lbs) weight rating



ALPS® Lock with Pyramid. Part Number and Spare Parts

Spare parts	Description	
S498-P	Main body with pyramid	
LDH - 694 A	Locking mechanism	
LDP - 695	Serrated pin	
LDA - 694	Hex nut adapter - Dummy	

Product Components





Lock with Flat Bottom

The ALPS® Lock Flat Bottom has shown proven reliability and smooth operation. This lightweight locking device made of glass reinforced nylon exhibits better wear characteristics than delrin or aluminum.

The convenient one-way gear allows the user to push in or wind out the limb into the prosthetic socket.

A unique feature of this lock is that if foreign material or clothing get caught in the gear, the housing can be removed, freeing the pin.



Key Benefits



Infinite pin adjustment reduces play



The housing can be removed, freeing the pin



For laminated or thermoplastic sockets



Lightweight (107gr)



No weight limit



ALPS® Lock Flat Bottom – Without Plate (S496-W) Part Number and Spare Parts

Spare parts	Description		
S496-W	Main body with flat bottom (without plate)		
LDH - 694 A	Locking mechanism		
LDP - 695	Serrated pin		
LDA - 694	Hex nut adapter – Dummy		



1 STAINLESS STEEL SERRATED PIN 3/8 x 2-5/8" LONG 1/4-20 US STANDARD THREADS PN# LDP-695



1 HEX NUT ADAPTER 1/4-20 US STANDARD THREADS TO 10mm PN# LDA-694



1 LOCKING MECHANISM (CLUTCH) PN# LDH-694A

ALPS® Lock Flat Bottom – With Plate (S496-T) Part Number and Spare Parts

Spare parts	Description		
S496-T	Main body with flat bottom (with plate)		
LDH - 694 A	Locking mechanism		
LDP - 695	Serrated pin		
LDA - 694	Hex nut adapter - Dummy		
UMP-496	Universal Mounting Plate		



1 STAINLESS STEEL SERRATED PIN 3/8 x 2-5/8" LONG 1/4-20 US STANDARD THREADS PN# LDP-695



1 HEX NUT ADAPTER 1/4-20 US STANDARD THREADS TO 10mm PN# LDA-694



1 LOCKING MECHANISM (CLUTCH) PN# LDH-694A



x1 UNIVERSAL MOUNTING PLATE PN# UMP-496



x6 #6 STAINLESS STEEL SCREWS SOLD WITH MOUNTING PLATE




x4 3/4" STAINLESS STEEL SOCKET HEAD CAP SCREWS SOLD WITH MOUNTING PLATE



ADDED SOLUTIONS





Prosthetic Socks

ALPS® prosthetic socks are light and soft, extremely comfortable with a perfect fit. They are made with technologically advanced fibers, highly effective in moisture management, that wick away perspiration so it can easily evaporate. This advanced fabric means that socks maintain their original quality for a long time, without getting misshapen or matting up.

Available with without distal hole (N), with distal hole (H), or with a reinforced distal hole (HR):









Choosing the best sock to suit your needs is easy with ALPS® socks' color-coded thickness indicator:

1 ply



3 ply (



5 ply



Key Benefits



Assists with volume fluctuation



Keeps its shape over time; minimal stretch and will not mat down



Technical fibers allow a cool, comfortable wear



Light and soft, extremely comfortable



SOCKS ARE AVAILABLE IN BUNDLES!


Contact your sales representative for details on all available bundle options.



Prosthetic Socks KC. Sizing Chart.

DIRECTIONS FOR SIZING

- 1. Measure from the distal end of the residual limb to the top of the prosthetic socket. The length of the correct sock should be approximately 2 inches more than this measurement.
- 2. Measure the "distal circumference" 2 inches up from the bottom of the distal end.



NARROW (KCN) SIZE	HOLE	LENGTH	1-PLY (1.5mm)	3-PLY (2.0mm)	5-PLY (2.5mm)
8"-9" No Hole Distal (N)		5"-7" (Extra Short)	KCN 08-1N	KCN 08-3N	KCN 08-5N
	No Hole	7"-9" (Short)	KCN 10-1N	KCN 10-3N	KCN 10-5N
	(N)	9"-15" (Regular)	KCN 14-1N	KCN 14-3N	KCN 14-5N
		15"-19" (Long)	KCN 18-1N	KCN 18-3N	KCN 18-5N
8"-9" Hole with Distal Reinforcement Circumference (HR)		5"-7" (Extra Short)	KCN 08-1HR	KCN 08-3HR	KCN 08-5HR
		7"-9" (Short)	KCN 10-1HR	KCN 10-3HR	KCN 10-5HR
		9"-15" (Regular)	KCN 14-1HR	KCN 14-3HR	KCN 14-5HR
		15"-19" (Long)	KCN 18-1HR	KCN 18-3HR	KCN 18-5HR
		5"-7" (Extra Short)	KCN 08-1H	KCN 08-3H	KCN 08-5H
	Hole without Reinforcement (H)	7"-9" (Short)	KCN 10-1H	KCN 10-3H	KCN 10-5H
		9"-15" (Regular)	KCN 14-1H	KCN 14-3H	KCN 14-5H
		15"-19" (Long)	KCN 18-1H	KCN 18-3H	KCN 18-5H

MEDIUM (KCM) SIZE	HOLE TYPE	LENGTH OF SOCKET	1-PLY (1.5mm)	3-PLY (2.0mm)	5-PLY (2.5mm)
9"-12" No Ho Distal (N)		5"-7" (Extra Short)	KCM 08-1N	KCM 08-3N	KCM 08-5N
	No Hole	7"-9" (Short)	KCM 10-1N	KCM 10-3N	KCM 10-5N
	(N)	9"-15" (Regular)	KCM 14-1N	KCM 14-3N	KCM 14-5N
		15"-19" (Long)	KCM 18-1N	KCM 18-3N	KCM 18-5N
Distal Reinforce		5"-7" (Extra Short)	KCM 08-1HR	KCM 08-3HR	KCM 08-5HR
	Hole with	7"-9" (Short)	KCM 10-1HR	KCM 10-3HR	KCM 10-5HR
	(HR)	9"-15" (Regular)	KCM 14-1HR	KCM 14-3HR	KCM 14-5HR
		15"-19" (Long)	KCM 18-1HR	KCM 18-3HR	KCM 18-5HR
Distal Reinfor		5"-7" (Extra Short)	KCM 08-1H	KCM 08-3H	KCM 08-5H
	Hole without Reinforcement (H)	7"-9" (Short)	KCM 10-1H	KCM 10-3H	KCM 10-5H
		9"-15" (Regular)	KCM 14-1H	KCM 14-3H	KCM 14-5H
		15"-19" (Long)	KCM 18-1H	KCM 18-3H	KCM 18-5H

WIDE (KCW) SIZE	HOLE TYPE	LENGTH OF SOCKET	1-PLY (1.5mm)	3-PLY (2.0mm)	5-PLY 2.5mm)
12"-14" Distal Circumference		5"-7" (Extra Short)	KCW 08-1N	KCW 08-3N	KCW 08-5N
	No Hole (N)	7"-9" (Short)	KCW 10-1N	KCW 10-3N	KCW 10-5N
		9"-15" (Regular)	KCW 14-1N	KCW 14-3N	KCW 14-5N
		15"-19" (Long)	KCW 18-1N	KCW 18-3N	KCW 18-5N
Distal Reinfor		5"-7" (Extra Short)	KCW 08-1HR	KCW 08-3HR	KCW 08-5HR
	Hole with	7"-9" (Short)	KCW 10-1HR	KCW 10-3HR	KCW 10-5HR
	Reinforcement (HR)	9"-15" (Regular)	KCW 14-1HR	KCW 14-3HR	KCW 14-5HR
		15"-19" (Long)	KCW 18-1HR	KCW 18-3HR	KCW 18-5HR
12"-14" Distal I Circumference		5"-7" (Extra Short)	KCW 08-1H	KCW 08-3H	KCW 08-5H
	Hole without Reinforcement (H)	7"-9" (Short)	KCW 10-1H	KCW 10-3H	KCW 10-5H
		9"-15" (Regular)	KCW 14-1H	KCW 14-3H	KCW 14-5H
		15"-19" (Long)	KCW 18-1H	KCW 18-3H	KCW 18-5H



Skin Reliever

The ALPS® Skin Reliever eliminates shear forces, friction, and prevents abrasions on the skin of the stump.

It also accommodates volume fluctuations over the life of the prosthesis.

There is no need to make a new socket, ALPS® Skin Reliever can be worn under silicone, gel, and pelite liners.



Key Benefits



Reduces shearing and abrasions



High elasticity



Replaces pelite and socks



Softens and soothes the skin



Contains an antioxidant for use on delicate skin




Skin Reliever ENCP. Sizing Chart.

Measure the circumference of the stump at 6 cm above the distal end.

MEASURED SIZE (CM)	ENCP 3/1.5MM TAPERED
10-19 cm	ENCP16
20-23 cm	ENCP20
24-27 cm	ENCP24
28-37 cm	ENCP28

Cross section





Antioxidant Cream

Patent No. 10,543,167 | 11,478,419

ALPS® Antioxidant Cream contains antioxidants that are effective scavengers of free radicals and hydroxyl groups. May assist with soothing pain, itchiness and irritated skin. Can be applied on stumps with sensitive skin as it stimulates skin regeneration process. Its moisturizing power improves skin elasticity.

Recommended for use with all liners.



Key Benefits



Soothes Pain, Itch, and Irritation



Stimulates Skin Regeneration



Improves Skin Elasticity



Great for Sensitive Skin

PRODUCT NO.	AMOUNT	QUANTITY
ANT-C50	50 ML	x 1
ANT-C50 BX12	50 ML	BOX 0F 12



Silicone Skin Lotion



ALPS® 100% Silicone Skin Lotion helps to soothe dry, sensitive skin. It is hypoallergenic, with no perfumes or dyes, and minimal residue. Can be used sparingly to relieve dry, cracked skin conditions.

Key Benefits



Hypoallergenic and skin friendly: formulated to soothe and hydrate sensitive skins without oily feeling



Simple use: apply a few drops of lotion on the skin or on palm of the hand and massage until completely absorbed



Concentrated: a few drops are enough to relieve dry and cracked skin conditions



Without soap or dyes: 100% free of soaps and dyes that could be aggressive

PRODUCT NO.	AMOUNT	QUANTITY
SL108	118 ML	x 1
SL108 BX12	118 ML	BOX 0F 12
SL108-32	944 ML	x 1



Prosthetic Lubricant

ALPS® Prosthetic Lubricant is used to help in donning gel and silicone prosthetic liners. It allows the amputee to don the liner or sheath easily while eliminating the sticky feel.



Key Benefits



Helps the placement of the liner in the socket



Easy application in 3 fast steps: Turn the liner inside-out - spray the lubricant on the external surface - unroll the liner over the stump



Eliminates the sticky feeling and is recommended for use with all liners



Contains 70% isopropyl alcohol for quick evaporation

PRODUCT NO.	AMOUNT	QUANTITY
PA604	15 OZ	x 1
PA604 BX12	15 OZ	BOX 0F 12



Prosthetic Cleanser



ALPS® Prosthetic Cleanser is a non-alkaline, soap-free cleansing agent designed to help remove dirt and excess perspiration from prosthetic liners and sleeves. Its unique formula allows it to work efficiently with hand washing or machine wash settings to effectively clean and restore items back to their original state.

Key Benefits



Developed to remove dirt and excess sweat from prosthetic liners and sleeves, either by hand or in the washing machine. It is recommended for use with all liners



Suitable for sensitive skin: thanks to its delicate formula it can also be applied to the areas of the skin that come into contact with the prosthesis



Without soap or dyes: 100% free of soaps and dyes that could be aggressive



Balanced ph: the non-alkaline silicone-based formula soothes the skin and does not damage the liner.

PRODUCT NO.	AMOUNT	QUANTITY
PD595	16 OZ	x 1
PD595 BX12	16 OZ	BOX 0F 12



Prosthetic Anti-Perspirant*

The first anti-perspirant truly designed for prosthetics. Created to help manage the continuing problem of perspiration, it has been formulated with twice the Aluminum Chlorohydrate found in over the counter anti-perspirants. Recommended for use with all liners.



Key Benefits



Manages issues of continuous perspiration



Made with twice the amount of Aluminum Chlorohydrate



Simple to use, fine mist sprayer with clear dust cap



Recommended for use with all ALPS gel liners

Product Code

PRODUCT NO.	AMOUNT	QUANTITY
ATP008	4 OZ	x 1
ATP008 BX 12	4 OZ	BOX 0F 12

*NOT AVAILABLE OUTSIDE OF USA



Prosthetic Ointment*



ALPS Prosthetic Ointment is a vitamin A and D formulation that soothes and protects the skin. It may be applied nightly after the prosthesis has been removed. Ideal for amputees with sensitive skin or irritation. May not be used while wearing ALPS Gel Products. Recommended for use with Silicone-CLEARPRO and Urethane.

Key Benefits



Made with Vitamin A & D



Soothes and protects skin



Ideal for sensitive or irritated skin

Product Code

PRODUCT NO.	AMOUNT	QUANTITY
P0840	4 OZ	x 1
P0840 BX12	4 OZ	BOX 0F 12

*NOT AVAILABLE OUTSIDE OF USA



PadSil

The Padsil was developed to increase patient comfort in combination with ALPS® liners. The more active user will appreciate the additional cushioning, softness, and durability. Users who experience volume changes will also appreciate the Padsil's ability to lift the limb back into the proper position. They can be placed at the bottom of the socket or inside a pin locking liner.



Key Benefits



Soft medical grade silicone



Helps maintain total contact



Keeps its shape over time



Easily modified



Blue dummy padsils are available for laminating



PadSil PDS. Sizing chart.

Measure the circumference of the stump at 6 cm above the distal end.

MEASURED SIZE (CM)	SILICONE DISTAL PAD	BLUE DUMMY FOR LAMINATION
14-19 cm	PDS-S	PDSB-S
20-23 cm	PDS-M	PDSB-M
24-27 cm	PDS-L	PDSB-L
28-31 cm	PDS-XL	PDSB-XL
32-35 cm	PDS-XXL	PDSB-XXL

Cross section









BioStep® EVO

Patent No. 7,410,503~ | 8,236,062~ | 8,574,314~ | 8,808,395~

Ready for the next adventure.

The most performing among the BioStep® Family. An all-around foot, with extra shock absorption and increased safety on all terrains, providing a softer connection to the ground.

Suggested for demanding dynamic users searching for the best. Ideal also for hiking and trekking, as well for recreational sport and high impact activities.

Natural, responsive, gently supportive.



Activity Level









K1

<2

3 🗸

Key Benefits



Lightweight



Physiological and energy efficient gait



Smooth rollover during the whole gait cycle



High dynamic response to support an active lifestyle



Possibility to do a wide range of activities with just one single foot



Ground compliance and stability on all terrains



Reduction of secondary issues thanks to a more balanced posture and symmetrical gait





Size Range: 22-30

Build Height: 17cm | 7in

Heel height: 1cm | .4in

Weight: 840 g | 30oz (size 26 with footshell)

Max user weight:

SIZE (cm)	К3
22	225lb
23-24	249lb
25-30	324lb

Product code example

Model	Size	Side	Rate
BSEV0	25	L/R	42

Add -D for dark footshell



Footshell and spectra sock are included

Key Features



NRG Technology®



Multilayer



Dynamic Calf and Shin



Split Toe



Range of Motion 26°



Axial Rotation 12°



Increased Shock Absorption



Waterproof

FIND YOUR PRODUCT CODE!





BioStep® PRO

Patent No. 7,410,503~ | 8,236,062~ | 8,574,314~ | 8,808,395~

Do that extra mile!

All benefits of the NRG Technology® with enhanced torsional movement.

Suggested for active outdoor patients who prefer a firm connection to the ground and the added benefits of torsion. Also for occasional and recreational sport activities.

Natural, comfortable, responsive.



Activity Level









K1

(2

/

Key Benefits



Lightweight



Physiological and energy efficient gait



Smooth rollover during the whole gait cycle



High dynamic response to support an active lifestyle



Possibility to do a wide range of activities with just one single foot



Ground compliance and stability on all terrains



Reduction of secondary issues thanks to a more balanced posture and symmetrical gait





Size Range: 22-30

Build Height: 15cm | 5.9in

Heel height: 1cm | .4in

Weight: 730g | 26oz (size 26 with footshell)

Max user weight:

SIZE (cm)	К3
22-24	223lb
25-30	403lb

Product code example

Model	Size	Side	Rate
BSPR0	25	L/R	42

Add -D for dark footshell



Footshell and spectra sock are included

Key Features



NRG Technology®



Multilayer



Dynamic Calf and Shin



Split Toe



Range of Motion 26°



Axial Rotation 12°



Increased Shock
Absorption



Waterproof

FIND YOUR PRODUCT CODE!





BioStep® CL

Patent No. 7,410,503~ | 8,236,062~ | 8,574,314~ | 8,808,395~

Your everyday companion.

Our Classic foot with all benefits of the NRG Technology®; it provides a firm connection to the ground.

Suggested for urban walkers and aggressive walkers searching for a strong, durable and comfortable dynamic response foot.

Natural and responsive.



Activity Level









K

Κ2

K3 🗸

K4

Key Benefits



Lightweight



Physiological and energy efficient gait



Smooth rollover during the whole gait cycle



High dynamic response to support an active lifestyle



Possibility to do a wide range of activities with just one single foot



Ground compliance and stability on all terrains



Reduction of secondary issues thanks to a more balanced posture and symmetrical gait





Size Range: 22-30

Build Height: 15cm | 5.9in

Heel height: 1cm | .4in

Weight: 740g | 26oz (size 26 with footshell)

Max user weight:

SIZE (cm)	К3
22-24	223lb
25-30	403lb

Product code example

Model	Size	Side	Rate
BSCL	25	L/R	42

Add -D for dark footshell



Footshell and spectra sock are included

Key Features



NRG Technology®



Multilayer



Dynamic Calf and Shin



Split Toe



Range of Motion 26°



Axial Rotation 12°



Increased Shock
Absorption



Waterproof

FIND YOUR PRODUCT CODE!





BioStep® XT

Patent No. 7,410,503~ | 8,236,062~ | 8,574,314~ | 8,808,395~

The relentless energy.

Maximum energy storage and return for the most demanding patients.

Suggested for high impact activities equiring torsional movement and for occasional sports activities.

Natural, comfortable, extremely responsive.



Activity Level









K

Κ2

K3 **✓**

Κ4

Key Benefits



Lightweight



Physiological and energy efficient gait



Smooth rollover during the whole gait cycle



High dynamic response to support an active lifestyle



Possibility to do a wide range of activities with just one single foot



Ground compliance and stability on all terrains



Reduction of secondary issues thanks to a more balanced posture and symmetrical gait





Size Range: 22-30

Build Height: 18.5cm | 7.3in

Heel height: 1cm | .4in

Weight: 840g | 30oz (size 26 with footshell)

Max user weight:

SIZE (cm)	K3
22-30	309lb

Product code example

Model	Size	Side	Rate
BSXT	25	L/R	42

Add -D for dark footshell



Footshell and spectra sock are included

Key Features



NRG Technology®



Multilayer



Dynamic Calf and Shin



Split Toe



Range of Motion 26°



Axial Rotation 12°



Increased Shock Absorption



Waterproof

FIND YOUR PRODUCT CODE!





BioStep® LP

Patent No. 8,574,314~

Low build height for a great freedom.

The ideal solution for active users with long residual limbs. It provides secure contact with the ground, as well as responsive action and ground compliance.

Suggested for medium and high mobility users with long residual limbs, who want a responsive, all-around foot.

Responsive with low height.



Activity Level









K1

Κ2

K3 **✓**

Key Benefits



Lightweight



Physiological and energy efficient gait



Smooth rollover during the whole gait cycle



High dynamic response to support an active lifestyle



Possibility to do a wide range of activities with just one single foot



Ground compliance and stability on all terrains



Reduction of secondary issues thanks to a more balanced posture and symmetrical gait





Size Range: 22-30

Build Height: 6.8cm | 2.7in

Heel height: 1cm | .4in

Weight: 600g | 21.2oz (size 26 with footshell)

Max user weight:

SIZE (cm)	К3
22-24	194lb
25-30	324lb

Product code example

Model	Size	Side	Rate
BSLP	25	L/R	Е

Add -D for dark footshell



Footshell and spectra sock are included

Key Features



Split Toe



Range of Motion 26°



Increased Shock Absorption



Waterproof

FIND YOUR PRODUCT CODE!





ProStride™

An intuitive microprocessor prosthetic knee that effortlessly adapts to the patient's gait speed, enhancing mobility and independence.

ProStride[™] knee is an advanced artificial joint designed to provide users with an easy-to-use, safe, comfortable, and energy-efficient walking experience.

ProStride $^{\text{TM}}$ includes microprocessor technology and a polycentric design, helping users easily initiate the swing phase while conserving energy, adapting to a wide range of walking conditions, and continuously adjusting to users' gait over time.



Activity Level









K

K2

K3 🗸

Key Features and Benefits



Four-bar pneumatic design.

Allows the instant center of rotation (ICR) trajectory to mimic that of a normal gait pattern.



Microprocessor technology.

Includes **three smart sensors** that automatically calibrate to the patient's gait pattern, making adjustments quick and accurate.



Smart programming.

Thanks to a **user-friendly app**, programming the ProStride can be completed within minutes with just a few simple inputs.



Easy to use.

The ProStride is designed to be user-friendly, enabling practitioners of all experience levels to install and program it efficiently, while users can easily manage it.



Customization options.


In the app, users have the opportunity to easily customize the knee's response to meet a specific activity need (biking, hiking, etc.) without affecting the initial programming settings.



Specifications

MATERIAL	Aluminum alloy (mainly), Stainless steel
SYSTEM WEIGHT	2.78lbs
MAX FLEXION ANGLE	142°
MAX BODY WEIGHT	275lbs
PROXIMAL CONNECTION	For connection with standard 4-hole adaptor
PYLON DIAMETER	1.18in
ACTIVITY LEVEL	K3
COMPATIBLE SOFTWARE	Mobile app compatible with Android Version 9+ and iOS Version 13+
BATTERY LIFE	Up to 3 days under normal usage
L CODES	L5840 and L5857

Dimensions



- A System Height
- **B** Alignment Reference Point
- **C** Proximal Build Height Up to the Alignment Reference Point
- **D** Distal Build Height Up to the Alignment Reference Point
- E Tube Insertion Length
- F Centro-Bolt to Alignment Reference

Note: All measurements are in mm





LINERS			
ITEM	DESCRIPTION	WARRANTY (SINGLE)	WARRANTY (PAIR)
AKDT/AKDTHD/ AKFR/AKFRHD	Extreme™ Liner	12 Months	24 Months
APDT-TP/APFR-TP	Anterior-Posterior Tapered Liner	12 Months	24 Months
ECDTHD/ECFRHD	Eco [™] Liner	12 Months	24 Months
ELDT/ELFR	EasyLiner [®]	12 Months	24 Months
ELPX	EasyLiner® Super Stretch	6 Months	
EZF	EasyFlex™ Liner	12 Months	24 Months
GPDT/GPFR	General Purpose Liner	12 Months	24 Months
PSI	Proximal Seal Interface Liner	6 Months	12 Months
SAFR/SADT/SBFR	Smart Seal™ Liner	6 Months	
SPDT/SPDTHD/SPFR/SPFRHD	Superior Performance Liner	12 Months	24 Months
VSDT	Liberty™ Liner	12 Months	24 Months
LSL/LSC	Silicone Pro Liner	12 Months	24 Months
PFGDT/PFGFR	Profit Liner	12 Months	24 Months
SSGSL/SSGSC	SoftSil® Liner	12 Months	24 Months
USHSL/USHSC/USLSL/USLSC/USUSL/USUSC	UltraSeal Silicone Liner	6 Months	
USHGL/USHGC/USLGL/USLGC/USUGL/USUGC	UltraSeal Gel Liner	6 Months	A



SLEEVES		
ITEM	DESCRIPTION	WARRANTY
SFB/SFBHD	Superior Suspension Sleeve	3 Months
SFK	Extreme Sleeve	3 Months
SFP/SFPHD	Superior Performance Sleeve	3 Months
SFR	Easy Sleeve	3 Months
SFX	Flex Sleeve	3 Months
SLE	Easy Sleeve Super Stretch	3 Months
VIVA	Valve Integrated Vacuum Activated	3 Months
RFS	Fabric Reinforced Straight Sleeve	3 Months
RFX®	Fabric Reinforced Flex Sleeve	3 Months
APS	Anterior Posterior Tapered Sleeve	3 Months

LOCKS		
ITEM	DESCRIPTION	WARRANTY
S498-P	ALPS® Lock with Pyramid	6 Months
S496-W	ALPS® Lock Flat Bottom	6 Months

EXTRA		
ITEM	DESCRIPTION	WARRANTY
ENCP	Skin Reliever	1 Month
KC	Prosthetic Socks	1 Month
PDS	Silicone Distal Pad	6 Months
VIP®	Vacuum Integrated Pump	6 Months

SKIN CARE		
ITEM	DESCRIPTION	WARRANTY
ANT-CEP50	Antioxidant Cream	90 Days
PA604	Prosthetic Lubricant	90 Days
PD595	Prosthetic Cleanser	90 Days
SL108	Skin Lotion	90 Days

BIOSTEP®	
ITEM	WARRANTY
BioStep® EVO, BioStep® PRO, BioStep® CL, BioStep® XT, BioStep® LP	3 Years
Foot Shell	6 Months
Spectra Sock	1 Month

KNEES	
ITEM	WARRANTY
ProStride™	3 Years

^{*}BioStep and ProStride do not have a 30 day no questions asked return policy.



CONTACT US

ALPS® WORLDWIDE OFFICES

ALPS® USA

Tel **727 528 8566** Tel **800 574 5426** Fax 727 528 8862

infous@easyliner.com www.easyliner.com

2895 42nd Avenue North St. Petersburg, FL 33714

ALPS® ITALY

Tel +39 051 86 06 42 info@alpsitaly.com www.alpsitaly.com

ALPS® EUROPE

Tel +420 377 223 127 info@easyliner.eu www.easyliner.eu

ALPS® UKRAINE

Tel +38(067) 889 80 22 alpsukraine@easyliner.eu www.alpsukraine.com

ALPS® CHINA


Tel 0731-55582976-801 china_salesrep@easyliner.com.cn www.easyliner.com.cn





For ALPS USA/Canada product inquiries, contact our sales team.

Locate your designated ALPS® sales representative by state.





CANDICE EVANS

West Territory Sales Manager

cevans@easyliner.com Phone: 727-342-8113 Mobile: 727-641-4977



DAN SCHEUMANN

Midwest Territory Sales Manager

dscheumann@easyliner.com Phone: 727-342-8111

Phone: 727-342-8111 Mobile: 727-744-9364



BRYAN RIZZO

Northeast Territory Sales Manager

brizzo@easyliner.com Phone: 727-342-8112 Mobile: 727-260-2693



NICK NELSON

Southeast Territory Sales Manager

nnelson@easyliner.com Phone: 727-342-8110 Mobile: 727-543-4558



Use of TM and (R) herein is limited to trademark claims and registrations in the United States.