



#### Nimm Adjust AFO

Nimm Adjust is neat and handy, providing great support that helps users with drop foot return to an active life. It lifts the foot immediately after toe-off, ensuring users can avoid obstacles such as carpet edges, thresholds, or even minor surface irregularities.

Neat and handy aid for drop foot patients



**Crafted from genuine leather:** soft on the skin and becoming more comfortable with daily use. Additional pads at the tibia enhance the orthosis' comfort.

**Easy to Attach**: The brace wraps securely around the ankle with Velcro straps. Its universal construction allows users to choose their preferred closure direction.

**Adjustable Lift**: The lift intensity can be adjusted in two ways: by repositioning the Velcro pad or relocating the hooks in the shoe.

**Automatic Foot Lift**: When walking, the elastic band automatically lifts the foot. Nimm Adjust is a discreet aid suitable for both children and adults.

### **NIMM ADJUST**



- Neat and handy aid for drop foot users
- Lifts the foot immediately after toe-off, preventing the user from worrying about tripping
- Made of leather that becomes softer with daily use
- Extra padding at tibia makes
  Nimm even more comfortable
- Universal design allows the user to choose whether to fasten the Velcro straps inward or outward
- Nimm Adjust is suitable for both children and adults





# Nimm Adjust AFO

Art.No.	Product	Size	L/R	
NA2-S	Nimm Adjust	Small/ barn 12-17 cm	Bilateral	
NA2-M	Nimm Adjust	Medium 15-20 cm	Bilateral	
NA2-L	Nimm Adjust	Large 20-27 cm	Bilateral	
NA2-XL	Nimm Adjust	X-Large >26 cm	Bilateral	





#### **Nimm Indoor**

Art.No.	Product	Size		
NI3-S	Nimm Indoor	Small/ barn	And Saint	
NI3-R	Nimm Indoor	Höger		9
NI3-L	Nimm Indoor	Vänster		

## **Nimm Adjust parts**

Art.No.	Product	Size	
NAE-S/M	Nimm Adjust Elastikgummiband	S/M 35 CM	
NAE-L/XL	Nimm Adjust Elastikgummiband	L/XL 40 CM	
			NIN
NSH	Nimm S-krokar NSH /pair		
NPH	Nimm plastkrokar /ea		43

