




Nimm Indoor

Many patients with foot drop sometimes want to be able to walk barefoot, especially indoors. Nimm Indoor makes it possible.

Nimm Indoor is a separate part to be used together with Nimm Regular or Nimm Adjust.

Nimm Indoor consists of a soft, durable leather cuff with two plastic hooks, the hooks are placed forward for greater lifting capacity. The inside has a friction pad that prevents the leather part from sliding and also adds comfort.



For details about the ankle part – see Nimm Regular and Nimm Adjust.



NIMM INDOOR



- Slim and easy aid for drop foot patients
- Made of leather that becomes softer with each day of use
- Nimm Indoor makes it possible to walk barefoot
- Nimm Indoor is suitable for both children and adults

Item.Nr.	Product	Size	 
NI3-S	Nimm Indoor	Small/ child	
NI3-R	Nimm Indoor	Right	
NI3-L	Nimm Indoor	Left	